



**YOUR PAIN RELIEF  
COMFORTABLY, AT YOUR HOME.**





## ABOUT

We are specialists in pain therapy and home rehabilitation.



## OUR MISSION

Improving patients' quality of life by reducing pain and increasing mobility.



## WHAT WE DO FOR YOU

We assist you during all your rehabilitation path: from the beginning until you will feel better. We also support you in the use of your device, providing simple tips.

### WE MANUFACTURE ELECTROMEDICAL DEVICES



**Safe** and easy to use at home



With **preset** therapeutic programs



Clinically approved and **certified** by bodies notified by the Ministry of Health, TÜV Rheinland Italia and M.T.I.C. Intercert



From the **design** and **development** of the products to the **customer care** service

**IMPORTANT:** We recommend that you always consult your doctor before and after treatment.

# MAGNETOTHERAPY

Magnetotherapy reduces pain and inflammation in bones, muscles and joints, due to traumatic and chronic pathologies. The technique uses the emission of pulsed electromagnetic fields (CEMP) that help to restore cell balance, repairing damaged cells.



## IDEAL FOR

- Osteoporosis
- Arthrosis
- Arthritis
- Inflammation
- Bone oedema
- Fractures
- Low back pain
- Algodistrophy
- Joint pain
- Ulcers
- Neuropathies


## USEFUL TIPS FOR A SUCCESSFUL MAGNETOTHERAPY TREATMENT

- ✓ Consult your doctor to identify the pathology
- ✓ Perform daily cycles of a minimum duration of 2-4 hours, depending on the pathology
- ✓ Duration of treatment cycles: 45-90 days (or greater)
- ✓ If necessary, repeat treatment throughout the year (i.e. in case of osteoporosis, arthrosis, arthritis)

# LAMAGNETO

LaMagnetto is the latest model of magnetotherapy that allows you to reduce pain and inflammation comfortably at home.

It features 46 pre-set therapeutic programs, updated according to the latest clinical trials, and 2 independent channels to carry out two different treatments at the same time.

 **Device designed and certified in Italy**

## SPECIFICATIONS

- ✓ 46 preset therapeutic programs
- ✓ 2 independent channels
- ✓ Colour display
- ✓ Intensity up to 150 gauss per channel
- ✓ Frequency from 1 to 100 Hz



**RECOMMENDED  
BY DOCTORS**

**INCLUDED WITH THE EQUIPMENT**



**46 THERAPEUTIC PROGRAMS  
UPDATED ACCORDING TO THE  
LATEST CLINICAL TRIALS**

**2 TREATMENTS  
SIMULTANEOUSLY**

**THE ONLY MAGNETOTHERAPY  
CERTIFIED FOR  
NEUROPATHY**

## PROGRAMS

<b>01</b> Osteoporosis	<b>17</b> Coxarthrosis	<b>33</b> Neuropathy
<b>02</b> Arthrosis	<b>18</b> Muscle atrophy	<b>34</b> Muscle strain
<b>03</b> Arthritis	<b>19</b> Muscle treatment	<b>35</b> Muscle cramp
<b>04</b> Cervical arthrosis	<b>20</b> Osteonecrosis	<b>36</b> Rhizarthrosis
<b>05</b> Joint pain	<b>21</b> Algodystrophy	<b>37</b> Impingement syndrome
<b>06</b> Neck pain	<b>22</b> Cartilage injury	<b>38</b> Carpal tunnel syndrome
<b>07</b> Chronic pain	<b>23</b> Ligament injury	<b>39</b> Titanium prosthesis
<b>08</b> Fractures	<b>24</b> Bone oedema	<b>40</b> Rotator cuff
<b>09</b> Epicondylitis	<b>25</b> Tendinitis	<b>41</b> Tarsal tunnel
<b>10</b> Epicondylitis (Golfer's elbow)	<b>26</b> Chondropathy	<b>42</b> Night 1 Hz
<b>11</b> Pseudoarthrosis	<b>27</b> Anti-inflammatory	<b>43</b> Night 15 Hz
<b>12</b> Low back pain	<b>28</b> Whiplash	<b>44</b> Night 50 Hz
<b>13</b> Lumbar-sacral pain	<b>29</b> Scarring	<b>45</b> Night 75 Hz
<b>14</b> Shoulder arthrosis	<b>30</b> Skin ulcers	<b>46</b> Automatic scanning (Autoscan)
<b>15</b> Knee arthrosis	<b>31</b> Discopathy	
<b>16</b> Scapulohumeral periarthritis	<b>32</b> Muscle pain	

## LAMAGNETO PRO

LaMagneto Pro is the ideal programmable magnetotherapy device for those who have a family member with difficulty in using electromedical devices.

 **Device designed and certified in Italy**

### SPECIFICATIONS

- ✓ 63 preset therapeutic programs  
+ 10 free customizable memories
- ✓ 2 independent channels
- ✓ Colour display
- ✓ Intensity up to 200 gauss per channel
- ✓ Frequency from 1 to 120 Hz
- ✓ **Counter of performed therapies**
- ✓ **Password protected therapies and settings**




INCLUDED WITH THE EQUIPMENT




## LAMAGNETO X

LaMagneto X is the ideal magnetotherapy device for night treatments all over the body.

 **Device designed and certified in Italy**

### SPECIFICATIONS

- ✓ 63 preset therapeutic programs  
+ 10 free customizable memories
- ✓ 2 independent channels
- ✓ Colour display
- ✓ Intensity up to 200 gauss per channel
- ✓ Frequency from 1 to 120 Hz
- ✓ Counter of performed therapies
- ✓ Password protected therapies and settings
- ✓ **Includes the Osteomat mattress for night or long lasting treatments all over the body**



INCLUDED WITH THE EQUIPMENT




## PROGRAMS

<b>01</b> Osteoporosis	<b>17</b> Coxarthrosis	<b>33</b> Neuropathy	<b>49</b> Treatment 5 Hz
<b>02</b> Arthrosis	<b>18</b> Muscle atrophy	<b>34</b> Muscle strain	<b>50</b> Treatment 8 Hz
<b>03</b> Arthritis	<b>19</b> Muscle treatment	<b>35</b> Muscle cramp	<b>51</b> Treatment 10 Hz
<b>04</b> Cervical arthrosis	<b>20</b> Osteonecrosis	<b>36</b> Rhizarthrosis	<b>52</b> Treatment 15 Hz
<b>05</b> Joint pain	<b>21</b> Algodystrophy	<b>37</b> Impingement syndrome	<b>53</b> Treatment 16 Hz
<b>06</b> Neck pain	<b>22</b> Cartilage injury	<b>38</b> Carpal tunnel syndrome	<b>54</b> Treatment 20 Hz
<b>07</b> Chronic pain	<b>23</b> Ligament injury	<b>39</b> Titanium prosthesis	<b>55</b> Treatment 30 Hz
<b>08</b> Fractures	<b>24</b> Bone oedema	<b>40</b> Rotator cuff	<b>56</b> Treatment 40 Hz
<b>09</b> Epicondylitis	<b>25</b> Tendinitis	<b>41</b> Tarsal tunnel	<b>57</b> Treatment 50 Hz
<b>10</b> Epicondylitis (Golfer's elbow)	<b>26</b> Chondropathy	<b>42</b> Night 1 Hz	<b>58</b> Treatment 60 Hz
<b>11</b> Pseudoarthrosis	<b>27</b> Anti-inflammatory	<b>43</b> Night 15 Hz	<b>59</b> Treatment 70 Hz
<b>12</b> Low back pain	<b>28</b> Whiplash	<b>44</b> Night 50 Hz	<b>60</b> Treatment 75 Hz
<b>13</b> Lumbar-sacral pain	<b>29</b> Scarring	<b>45</b> Night 75 Hz	<b>61</b> Treatment 80 Hz
<b>14</b> Shoulder arthrosis	<b>30</b> Skin ulcers	<b>46</b> Automatic scanning (Autoscan)	<b>62</b> Treatment 90 Hz
<b>15</b> Knee arthrosis	<b>31</b> Discopathy	<b>47</b> Treatment 1 Hz	<b>63</b> Treatment 100 Hz
<b>16</b> Scapulohumeral periarthritis	<b>32</b> Muscle pain	<b>48</b> Treatment 3 Hz	

# MAG2000

Mag2000 is the world's best-selling magnetotherapy device.



 Device designed and certified in Italy

## SPECIFICATIONS

- ✓ 35 preset therapeutic programs
- ✓ 2 channels
- ✓ Intensity up to 150 gauss per channel
- ✓ Frequency from 1 to 100 Hz

INCLUDED WITH THE EQUIPMENT



## PROGRAMS


01 Osteoporosis	10 Epicondylitis (Golfer's elbow)	19 Muscle contractures	28 Treatment 40 Hz
02 Arthrosis	11 Intercostal bruising	20 Osteonecrosis	29 Treatment 50 Hz
03 Arthritis	12 Low back pain	21 Treatment 1 Hz	30 Treatment 60 Hz
04 Cervical arthrosis	13 Lumbar-sacral pain	22 Treatment 3 Hz	31 Treatment 70 Hz
05 Joint pain	14 Shoulder arthrosis	23 5 Hz treatment	32 Treatment 80 Hz
06 Neck pain	15 Knee arthrosis	24 Treatment 10 Hz	33 Treatment 90 Hz
07 Sprains	16 Scapulohumeral periarthritis	25 Treatment 15 Hz	34 Treatment 100 Hz
08 Fractures	17 Coxarthrosis	26 Treatment 20 Hz	35 Automatic scanning (Autoscan)
09 Epicondylitis	18 Muscle atrophy	27 Treatment 30 Hz	



# MAG2000 PLUS

Mag2000 Plus is the ideal magnetotherapy device for high-intensity treatments or in presence of plaster.



 **Device designed and certified in Italy**

## SPECIFICATIONS

- ✓ 35 preset therapeutic programs  
+ 10 free customizable memories
- ✓ 2 channels
- ✓ Intensity up to 200 gauss per channel
- ✓ Frequency from 1 to 120 Hz

INCLUDED WITH THE EQUIPMENT



## PROGRAMS

01 Osteoporosis	10 Epicondylitis (Golfer's elbow)	19 Muscle contractures	28 Treatment 40 Hz
02 Arthrosis	11 Intercostal bruising	20 Osteonecrosis	29 Treatment 50 Hz
03 Arthritis	12 Low back pain	21 Treatment 1 Hz	30 Treatment 60 Hz
04 Cervical arthrosis	13 Lumbar-sacral pain	22 Treatment 3 Hz	31 Treatment 70 Hz
05 Joint pain	14 Shoulder arthrosis	23 5 Hz treatment	32 Treatment 80 Hz
06 Neck pain	15 Knee arthrosis	24 Treatment 10 Hz	33 Treatment 90 Hz
07 Sprains	16 Scapulohumeral periarthritis	25 Treatment 15 Hz	34 Treatment 100 Hz
08 Fractures	17 Coxarthrosis	26 Treatment 20 Hz	35 Automatic scanning (Autoscan)
09 Epicondylitis	18 Muscle atrophy	27 Treatment 30 Hz	

# MAG700

Mag700 is the one-channel magnetotherapy device capable of acting on some osteoarticular pathologies.

 **Device designed and certified in Italy**



## SPECIFICATIONS

- ✓ 11 programs
- ✓ 1 channel
- ✓ Intensity up to 70 gauss
- ✓ Frequency from 5 to 100 Hz

INCLUDED WITH THE EQUIPMENT




## PROGRAMS

01 Osteoarthritis	04 Fractures	07 Anti-inflammatory	10 Treatment 10 Hz
02 Arthritis	05 Joint pain	08 Automatic scanning (Autoscan)	11 Treatment 20 Hz
03 Osteoporosis	06 Muscle pain	09 5 Hz treatment	

## OSTEOMAT

The Osteomat mattress is the most comfortable and practical solution for night treatments all over the body.

 **Device designed and certified in Italy**

### COMPATIBLE WITH

LaMagneto, LaMagneto Pro, LaMagneto X, Mag2000, Mag2000 Plus



**8 SOLENOIDS FOR A TRUE “FULL BODY” TREATMENT**

**SOFT AND FLUFFY MATERIAL WITH ULTRA-THIN SOLENOIDS**

**LIGHTWEIGHT AND FOLDABLE TO BE TRANSPORTED WHEREVER YOU WANT**

## TAP2000

40x40cm mat with 2 cables of 3 solenoids each (in total 6 solenoids).



### COMPATIBLE WITH

Mag2000, Mag2000 Plus



## CIL-MAG

Professional cylindrical solenoid (diameter 25cm), lightweight and transportable for limb treatment.



### COMPATIBLE WITH

Mag2000, Mag2000 Plus



## HP2S2000

Pair of professional solenoids for high intensity therapy (solenoids diameter: 11 cm). Suitable for treatments up to 150-200 gauss (minimizes solenoids heating) and for applications over plaster.

### COMPATIBLE WITH

LaMagneto, LaMagneto Pro, LaMagneto X, Mag2000, Mag2000 Plus



## F3S2000

Therapeutic elastic band for low frequency magnetotherapy with 3 solenoids cable.

### COMPATIBLE WITH

LaMagneto, LaMagneto Pro, LaMagneto X, Mag2000, Mag2000 Plus



## F2S700

Therapeutic elastic band for low frequency magnetotherapy with 2 solenoids cable.


### COMPATIBLE WITH

Mag700



## ORTHOMAG

Orthomag is the first low-frequency magnetotherapy device that can be used anywhere, in total freedom of movement.

 Device designed and certified in Italy

### SPECIFICATIONS

- ✓ 2 programs (50-75 Hz)
- ✓ Rechargeable Lithium Battery
- ✓ Up to 25 gauss per program



WORKS BOTH WIRED AND WIRELESS

GUARANTEES TOTAL FREEDOM OF MOVEMENT

INCLUDED WITH THE EQUIPMENT



FLEXIBLE  
APPLICATOR



For more information visit [www.orthomag.eu](http://www.orthomag.eu)

## MAG1000

Mag1000 is the high-frequency magnetotherapy device ideal for the treatment of soft tissue diseases such as muscles and tendons.

 **Device designed and certified in Italy**

### SPECIFICATIONS

- ✓ 20 preset therapeutic programs
- ✓ 2 channels
- ✓ Intensity adjustment on 3 levels
- ✓ Frequency from 8 to 640 Hz



**WORKS ON BATTERY**

**INCLUDED WITH THE EQUIPMENT**



### PROGRAMS

<b>01</b> Arthritis/ Osteoarthritis	<b>06</b> Osteoporosis	<b>11</b> Bruises	<b>16</b> Epicondylitis (Golfer's elbow)
<b>02</b> Rheumatism	<b>07</b> Fractures	<b>12</b> Myalgia	<b>17</b> Low back pain
<b>03</b> Pains/Joint Prosthesis	<b>08</b> Sprains	<b>13</b> Tendinitis	<b>18</b> Muscle contractures
<b>04</b> Herniated discs	<b>09</b> Dislocations	<b>14</b> Hematomas	<b>19</b> Average application 2h
<b>05</b> Muscle pain	<b>10</b> Traumas	<b>15</b> Epicondylitis	<b>20</b> 4h long application

### ACCESSORIES

#### FTI1000

Therapeutic elasticband for HF magnetotherapy with high frequency.



#### TAP-HF

40x40cm mat for high frequency magnetotherapy.



### HOW DOES MAGNETOTHERAPY WORK AGAINST PAIN?

The effectiveness of magnetotherapy is based on the pulsation of the electromagnetic field that is activated and deactivated alternately in rapid sequence. In this way, the cell is stimulated by the sequence of electromagnetic fields that give energy to the cell, restoring the electrical balance. Biological modifications induced by magnetic fields on biological membranes, cells and vessels have therapeutic effects such as: anti-inflammatory, anti-oedema and stimulate tissue repair. The relaxation of the muscles obtained with magnetotherapy promotes the improvement of microcirculation. Smooth fibers relax, with consequent vasodilation effect. The release of endorphins contributes to raise the pain threshold, increase the local blood flow and get an anti-inflammatory effect.

### ARE MAGNETOTHERAPY DEVICES EFFECTIVE?

Magnetotherapy has been a well-established therapy with universally recognised effectiveness for decades. There are thousands of clinical studies in the literature that report the efficacy of magnetotherapy in the treatment of pathologies of the musculoskeletal system.

### DOES MAGNETOTHERAPY HAVE CONTRAINDICATIONS?

Magnetotherapy is one of the physical therapies with the least number of contraindications as is also evidenced by the numerous clinical studies carried out on pulsed electromagnetic fields.

Some basic warnings indicated in the user manuals concerning: pregnancy, tuberculosis, juvenile diabetes, viral diseases in the acute phase, mycosis, heart disease, tumours, severe arrhythmias, patients with pacemakers, children, patients with magnetizable prostheses, acute infections, and epileptic patients. (Unless otherwise prescribed by your doctor).

### HOW LONG CAN YOU DO MAGNETOTHERAPY?

It is essential to consult your doctor first to properly diagnose the condition you are suffering from.

In general, magnetotherapy requires medium-long treatment times. It is recommended to perform cycles of 45-90 days (or higher) with a minimum duration of 2-4 hours a day, unless otherwise indicated by the doctor. If necessary, it is possible to repeat the treatment several times during the year. For example in the case of arthrosis, arthritis, osteoporosis and joint pain. There are no risks of overdose. In fact, as in the case of treatments for osteoporosis, it is possible to do night therapy even for 8 consecutive hours.

### WHAT IS THE DIFFERENCE BETWEEN LOW FREQUENCY AND HIGH FREQUENCY MAGNETOTHERAPY?

Low-frequency and high-intensity magnetotherapy consists of a device that emits pulsed electromagnetic fields (CEMP) transferred to the patient via solenoids.

CEMP are particularly suitable for the treatment of bone tissue pathologies, but important evidence is also found in the treatment of muscle-ligament problems.

High frequency and low intensity magnetotherapy, on the other hand, is based on the emission of short pulse radio waves with adjustable frequency between 18 MHz and 900 MHz. Clinical studies demonstrate its effectiveness for all soft tissue diseases, while maintaining important results also in the treatment of bone diseases.

# PRESSOTHERAPY

Pressotherapy is a medical device useful for the promotion of venous and lymphatic return of the lower and upper limbs.

It stimulates the circulatory system (arterial, venous and lymphatic) and drains body fluids, promoting the reactivation of venous circulation and the elimination of toxins.

## OPERATING PRINCIPLE

Pressotherapy supports the natural movement of flows involving the venous and lymphatic system.

It exerts controlled compression onto limbs and abdomen with a distal-proximal movement, that is, from the periphery to the center of the body.

Therapy parameters are set and activated by the airflow device at the suggested pressure, inflating and deflating intermittently and sequentially the applicator chambers.

Illness	Pressure (mmHg)	Time (minutes)
Oedema	50	50
Lymphedema	50	30
Venous ulcers	50	60
Venous insufficiency	50	60
Muscle recovery	80	20





## I-PRESS

I-Press is the new Made in Italy pressotherapy device, ideal for the treatment of circulatory system conditions. Thanks to its pumping action, it promotes venous return, increasing blood flow of tissues and their consequent proper physiological renewal.

**MADE IN ITALY**

### SPECIFICATIONS

- ✓ Nominal pressure: max 200 mmHg ( $\pm 20\%$ )
- ✓ Therapy time settable: from 0 to 30 minutes
- ✓ Maximum consumption: 100 mA
- ✓ Power supply: 230 VAC, 50 Hz
- ✓ Weight: 2.0 kg



**RECOMMENDED  
BY DOCTORS**



Accessories size guide from  
page 22 to page 24



1st clinical study in the world  
on home pressotherapy

## I-PRESS LEG1

- ✓ Machine body
- ✓ 1 leg (right or left)
- ✓ 4-termination tube
- ✓ 1 insole pressor



I-PRESS LEG1  
SIZE M



I-PRESS LEG1  
SIZE L



## I-PRESS LEG2

- ✓ Machine body
- ✓ 1 right leg
- ✓ 1 left leg
- ✓ 8-termination tube
- ✓ 2 insole pressor



I-PRESS LEG2  
SIZE M



I-PRESS LEG2  
SIZE L



## I-PRESS ARM1

- ✓ Machine body
- ✓ 1 armband (right or left)
- ✓ 4-termination tube



I-PRESS ARM1  
SIZE M



I-PRESS ARM1  
SIZE L



\* Supplied with armband + extender

All 3 configurations can be further customized with the entire range of pressotherapy accessories

All pictures shown are for illustration purpose only. They may vary in some physical aspects.

## POWER Q-1000 PREMIUM

Power Q-1000 Premium is a pressotherapy device with 6 modes of inflation, ideal for treating lymphedema, early cellulite, post-mastectomy arm drainage, arthritis rheumatoid, hematomas following surgical interventions and for postoperative prophylaxis in venous thrombosis interventions.



### SPECIFICATIONS

- ✓ Adjustable pressure: 0 to 240 mmHg ( $\pm 20\%$ )
- ✓ Adjustable pressure for each air chamber
- ✓ 9 levels of airflow speed
- ✓ 6 programs
- ✓ Therapy time settable: from 5 to 99 minutes
- ✓ Maximum consumption: 53 W
- ✓ Power supply: 220-240 VAC, 50-60 Hz
- ✓ Current: 0.2 TO
- ✓ Weight: 2.0 kg



Accessories size guide from  
page 22 to page 24

## TOTAL VERSION

- ✓ Machine body
- ✓ 2 4-chamber leggings
- ✓ 2 insoles
- ✓ 4 compression chambers abdominal band
- ✓ 4 compression chambers armband
- ✓ Adaptor for simultaneous leg and abdomen therapy
- ✓ Fitting hoses



POWER Q-1000  
PREMIUM  
SIZE M



POWER Q-1000  
PREMIUM  
SIZE L



## LEG2-ABD VERSION

- ✓ Machine body
- ✓ 2 4-chamber leggings
- ✓ 2 insoles
- ✓ 4 compression chambers abdominal band
- ✓ Adaptor for simultaneous leg and abdomen therapy
- ✓ Fitting hoses



POWER Q-1000  
PREMIUM  
SIZE M



POWER Q-1000  
PREMIUM  
SIZE L



## LEG2 VERSION

- ✓ Machine body
- ✓ 2 4-chamber leggings
- ✓ 2 insoles
- ✓ Fitting hose



POWER Q-1000  
PREMIUM  
SIZE M



POWER Q-1000  
PREMIUM  
SIZE L



## ARM1 VERSION

- ✓ Machine body
- ✓ 4 chamber armband
- ✓ Fitting hose



POWER Q-1000  
PREMIUM  
SIZE M



POWER Q-1000  
PREMIUM  
SIZE L



\* Supplied with armband + extender

All pictures shown are for illustration purpose only. They may vary in some physical aspects.

## LEGGINGS

4-chamber leggings

### LEGGINGS

Right & Left


THIGH CIRCUMFERENCE

ANKLE CIRCUMFERENCE

LEG LENGTH

(measurement taken from inside thigh to heel)

SIZE **(M)**




Up to 60 cm

Up to 33 cm

Up to 76 cm

SIZE **(L)**





Up to 72 cm

Up to 46 cm

Up to 84 cm

### EXTENDERS

Pair of leg extenders



### TNT PROTECTIONS

Pair of hygienic protections made of TNT for leggings



All pictures shown are for illustration purpose only. They may vary in some physical aspects.

ARBANDS

4-chamber armband

ARBAND


Right & Left

ARMPIT CIRCUMFERENCE

WRIST CIRCUMFERENCE

ARM LENGTH

(measure taken from the armpit to the wrist)



Up to 46 cm

Up to 32 cm

Up to 87 cm

EXTENDERS

Pair of armband extenders

ARMPIT CIRCUMFERENCE


(armband + extender)

WRIST CIRCUMFERENCE

(armband + extender)

ARM LENGTH

(measure taken from the armpit to the wrist)



Up to 59 cm

Up to 37 cm

Up to 87 cm


All pictures shown are for illustration purpose only. They may vary in some physical aspects.

## ABDOMINAL BAND

Abdominal band with 4 chambers


WAIST CIRCUMFERENCE

SIZE **(M)**



Up to 116 cm

SIZE **(L)**



From 117 to 132 cm

## ACCESSORIES

### DUAL CONNECTOR

Connector with 4 tubes and 8 terminations for simultaneous therapy (2 leggings or 2 armbands)



### SINGLE CONNECTOR

4 tube connector for armband, legging or abdominal band



### ADAPTER FOR DOUBLE THERAPY

Adaptor for simultaneous therapy (2 leggings and abdominal belt).



All pictures shown are for illustration purpose only. They may vary in some physical aspects.

# ELECTROTHERAPY

Electrotherapy consists in the transmission of micro electrical impulses to the human body for therapeutic purposes. Such pulses are characterized by precise parameters such as:


**Frequency:** the number of impulses per time unit, measured in Hertz (Hz)

**Pulse width:** the duration of the single pulse, measured in microseconds

**Intensity:** the amount of current induced in the body, measured in milliamperes (mA)

By modulating these parameters, several **therapeutic benefits can be obtained:**


- ✓ Muscle recovery
- ✓ Pain relief
- ✓ Treatment of incontinence
- ✓ Treatment of denervated muscle
- ✓ Fluid drainage





# I-TECH PHYSIO

I-Tech Physio is the easy-to-use and fully programmable electrostimulator for pain therapy and muscle rehabilitation.

 **Device designed and certified in Italy**

## SPECIFICATIONS

- ✓ 2 independent channels
- ✓ Compensated biphasic square wave
- ✓ Monophasic square wave
- ✓ Interferential wave
- ✓ Kotz
- ✓ Maximum Intensity: 100mA
- ✓ Ni-Mh rechargeable internal battery with LL technology (long-lasting) low self-discharge

**INCLUDES TREATMENTS FOR DENERVATED,  
IONTOPHORESIS AND INCONTINENCE**



**RECOMMENDED  
BY DOCTORS**


**INCLUDED WITH THE EQUIPMENT**



TENS	08 Tens Burst	03 Resistance of lower limbs	04 Tone-up of lower limbs
01 Tens rapid	09 Atrophy prevention	04 Strength resistance of upper limbs and trunk	05 Definition of upper limbs and trunk
02 Tens endorphinic	10 Atrophy	05 Strength resistance of lower limbs	06 Definition of lower limbs
03 Tens at maximum values	11 Denervated square 50 ms	06 Strength base upper limbs and trunk	07 Modelling
04 Anti-inflammatory	12 Denervated square 100 ms	07 Strength base lower limbs	08 Microlifting
05 Neck pain/Myotensive headache	13 Denervated square 150 ms	08 Fast strength upper limbs and trunk	09 Lipolysis abdomen
06 Low back pain/Sciatica	14 Denervated square 200 ms	09 Fast strength lower limbs	10 Lipolysis thighs
07 Distortions/Bruises	15 Denervated square 250 ms	10 Explosive force upper limbs and trunk	11 Lipolysis buttocks and hips
08 Vascularisation	16 Denervated triangular 50 ms	11 Explosive force lower limbs	12 Lipolysis arms
09 Muscle tension relieving	17 Triangular denervated 100 ms	12 Deep capillarisation	13 Tissue elasticity
10 Hand and wrist pain	18 Triangular denervated 150 ms	13 Muscle recovery	14 Capillarisation
11 Sole stimulation	19 Triangular denervated 200 ms	14 Agonist-antagonist	15 Heaviness in the legs
12 Epicondylitis	20 Denervated triangular 250 ms	15 Sequential tonic contractions upper limbs and trunk	INCONTINENCE
13 Epicondylitis (Golfer's elbow)	21 Interferential	16 Sequential tonic contractions lower limbs	01 Stress incontinence*
14 Periarthritis	22 Tens with amplitude modulation	17 Sequential phasic contractions upper limbs and trunk	02 Stress incontinence*
REHA	23 Tens alternating	18 Sequential phasic contractions lower limbs	03 Stress incontinence*
01 Iontophoresis L	24 Tens with frequency modulation	19 Muscle tension relieving	04 Urgency incontinence*
02 Iontophoresis M	25 Tens with frequency modulation	20 Deep massage	05 Urgency incontinence*
03 Iontophoresis H	26 Tens with amplitude modulation	21 EMS	06 Urgency incontinence*
04 Microcurrent	27 Kotz	BEAUTY	07 Mixed incontinence*
05 Haematoma	NEMS	01 Firming upper limbs and trunk	08 Mixed incontinence*
06 Oedema	01 Heating	02 Firming lower limbs	09 Mixed incontinence*
07 Sequential Tens	02 Resistance upper limbs and trunk	03 Toning upper limbs and trunk	*Probes not included


# MIO-CARE PRO

Mio-Care Pro is the ideal electrotherapy device for pain therapy and muscle stimulation.

 **Device designed and certified in Italy**

## SPECIFICATIONS

- 2 independent channels
- Compensated biphasic square wave
- Belt hook clip
- Maximum Intensity: 100mA
- Ni-Mh rechargeable internal battery with LL technology  
Low self-discharge (long-lasting)



INCLUDED WITH THE EQUIPMENT




## PROGRAMS

TENS	02	Resistance upper limbs and trunk	19	Muscle tension relieving	14	Capillarisation	
01	Tens rapid	03	Resistance of lower limbs	20	Deep massage	15	Heaviness in the legs
02	Tens endorphinic	04	Strength resistance of upper limbs and trunk	21	EMS	FREE MEMORIES	
03	Tens at maximum values	05	Strength resistance of lower limbs	BEAUTY		01	TENS Mem1
04	Anti-inflammatory	06	Strength base upper limbs and trunk	01	Firming upper limbs and trunk	02	TENS Mem2
05	Cervicalgia/Myotensive headache	07	Strength base lower limbs	02	Firming lower limbs	03	TENS Mem3
06	Low back pain/Sciatica	08	Fast strength upper limbs and trunk	03	Toning upper limbs and trunk	04	TENS Mem4
07	Distortions/Bruises	09	Fast strength lower limbs	04	Tone-up of lower limbs	05	TENS Mem5
08	Vascularisation	10	Explosive force upper limbs and trunk	05	Definition of upper limbs and trunk	06	NEMS Mem1
09	Muscle tension relieving	11	Explosive force lower limbs	06	Definition of lower limbs	07	NEMS Mem2
10	Hand and wrist pain	12	Deep capillarisation	07	Modelling	08	NEMS Mem3
11	Sole stimulation	13	Muscle recovery	08	Microlifting	09	NEMS Mem4
12	Epicondylitis	14	Agonist-antagonist	09	Lipolysis abdomen	10	NEMS Mem5
13	Epicondylitis (Golfer's elbow)	15	Sequential tonic contractions upper limbs and trunk	10	Lipolysis thighs	11	NEMS Alternating
14	Periarthritis	16	Sequential tonic contractions lower limbs	11	Lipolysis buttocks and hips	12	NEMS Alternating
NEMS	17	Sequential phasic contractions upper limbs and trunk	12	Lipolysis arms	13	Battery test	
01	Heating	18	Sequential phasic contractions lower limbs	13	Tissue elasticity		

## MIO-CARE FITNESS

Mio-Care Fitness is the ideal electrostimulator for the treatment of pain, able to meet the needs of the athletes and of those who want to improve their physical condition.

 **Device designed and certified in Italy**

### SPECIFICATIONS

- ✓ 2 independent channels
- ✓ Compensated biphasic square wave
- ✓ Belt hook clip
- ✓ Maximum Intensity: 100mA
- ✓ Ni-Mh rechargeable internal battery with LL technology  
Low self-discharge (long-lasting)




INCLUDED WITH THE EQUIPMENT



## MIO-CARE TENS

Mio-Care Tens is the best selling electrotherapy device for the treatment of pain.

 **Device designed and certified in Italy**

### SPECIFICATIONS

- ✓ 2 independent channels
- ✓ Compensated biphasic square wave
- ✓ Belt hook clip
- ✓ Maximum Intensity: 100mA
- ✓ Ni-Mh rechargeable internal battery with LL technology  
Low self-discharge (long-lasting)



INCLUDED WITH THE EQUIPMENT



## PROGRAMS


TENS (Mio-Care Tens and Mio-Care Fitness)			
01	Conventional Tens (Quick)	15	Epicondylitis
02	Endorphin Tens (delayed)	16	Epicondylitis (Golfer's elbow)
03	Tens at maximum values	17	Periarthritis
04	Anti-inflammatory	18	Microcurrent
05	Cervicalgia/Myotensive headache	19	Stress incontinence*
06	Low back pain/Sciatica	20	Urgency incontinence*
		FITNESS (Mio-Care Fitness)	
07	Distortions/Bruises	01	Heating
08	Vascularisation	02	Resistance upper limbs and trunk
09	Muscle tension relieving	03	Resistance of lower limbs
10	Hematomas	04	Strength resistance of upper limbs and trunk
11	Atrophy prevention	05	Strength resistance of lower limbs
12	Atrophy (trophism recovery)	06	Strength base upper limbs and trunk
13	Hand and wrist pain	07	Strength base lower limbs
14	Sole stimulation	08	Fast strength upper limbs and trunk
		BEAUTY (Mio-Care Fitness)	
09	Fast strength lower limbs	03	Toning upper limbs and trunk
10	Explosive force upper limbs and trunk	04	Tone-up of lower limbs
11	Explosive force lower limbs	05	Definition of upper limbs and trunk
12	Deep capillarisation	06	Definition of lower limbs
13	Muscle recovery	07	Modelling
14	Agonist-antagonist	08	Microlifting
15	Sequential tonic contractions upper limbs and trunk	09	Lipolysis abdomen
16	Sequential tonic contractions lower limbs	10	Lipolysis thighs
17	Sequential phasic contractions upper limbs and trunk	11	Lipolysis buttocks and hips
18	Sequential phasic contractions lower limbs	12	Lipolysis arms
19	Muscle tension relieving	13	Tissue elasticity
20	Deep massage	14	Capillarisation
		15	Heaviness in the legs

\*Probes not included

## MIO-IONOTENS

Mio-Ionotens is the electrostimulator dedicated to iontophoresis, a technique that uses direct current to introduce medicines in the painful area or contracture. In this way, it avoids the intake of oral medications.



 **Device designed and certified in Italy**

INCLUDED WITH THE EQUIPMENT

### PROGRAMS


- ✓ 3 IONTOPHORESIS programs
- ✓ 7 REHABILITATION programs
- ✓ 14 TENS programs
- ✓ 13 free TENS and NEMS memories



## MIO-PERISTIM

Mio-Peristim is the electrostimulation designed to treat problems related to male and female incontinence.

It stimulates the wall of the pelvic muscles, responsible for controlling the bladder, toning and strengthening it. Mio-Peristim increases muscle control and decreases bladder irritability.

 **Device designed and certified in Italy**

INCLUDED WITH THE EQUIPMENT

### PROGRAMS

- ✓ 9 programs INCONTINENCE
- ✓ 14 TENS programs
- ✓ 13 free TENS and NEMS memories



Vaginal probe



Anal probe



All pictures shown are for illustration purpose only. They may vary in some physical aspects.

# WHICH T-ONE TO CHOOSE?

T-ONE IS OUR RANGE OF 4-CHANNEL ELECTROSTIMULATORS.

T-ONE COACH



T-ONE MEDI SPORT



T-ONE MEDI PRO



T-ONE REHAB



IDEAL FOR	Amateur athletes	Habitual athletes	Sports medicine doctors	The rehabilitation specialists
TENS Programs	16	16	16	16
NEMS Programs	14	20	20	20
BEAUTY Programs	18	18	18	18
Programs for two users (2+2)	10	22	22	22
REHA Programs		10	16	29
URO Programs			9	9
Free memories TENS/NEMS		12	12	12



## T-ONE COACH

T-One Coach is the 4-channel electrostimulator designed for the prevention of trauma and muscle strengthening of amateur athletes.

 **Device designed and certified in Italy**

### SPECIFICATIONS

- ✓ 4 independent channels (up to 8 pluggable electrode)
- ✓ Constant current generator
- ✓ Compensated biphasic square wave
- ✓ Two-user function (2+2)
- ✓ Operation with 4.8 Volt rechargeable battery pack
- ✓ Adjustable intensity from 1 to 120 mA
- ✓ Wave frequency from 1 to 150 Hz
- ✓ Pulse width 20 to 300  $\mu$ s
- ✓ Preset time for each program



**INCLUDED WITH THE EQUIPMENT**




### PROGRAMS

- ✓ 16 TENS programs
- ✓ 14 NEMS programs
- ✓ 18 BEAUTY programs
- ✓ 10 double user programs (2+2)

For more information about the device and its programs: [www.t-onetherapy.com](http://www.t-onetherapy.com)

## T-ONE MEDI SPORT

T-One Medi Sport is the 4-channel electrostimulator dedicated to regular athletes who want to improve their performance and carry out rehabilitation independently.

 **Device designed and certified in Italy**

### SPECIFICATIONS

- ✓ 4 independent channels (up to 8 pluggable electrode)
- ✓ Constant current generator
- ✓ Compensated biphasic square wave
- ✓ Monophasic wave (only for iontophoresis programs)
- ✓ Two-user function (2+2)
- ✓ Operation with 4.8 Volt rechargeable battery pack
- ✓ Adjustable intensity from 1 to 120 mA
- ✓ Wave frequency from 1 to 1200 Hz
- ✓ Pulse width from 20 to 450  $\mu$ s
- ✓ Time from 1 to 90 minutes



**INCLUDED WITH THE EQUIPMENT**




### PROGRAMS

- |                      |                                 |
|----------------------|---------------------------------|
| ✓ 10 REHA programs   | ✓ 22 double user programs (2+2) |
| ✓ 16 TENS programs   | ✓ 12 TENS/NEMS free memories    |
| ✓ 20 NEMS programs   |                                 |
| ✓ 18 BEAUTY programs |                                 |

For more information about the device and its programs: [www.t-onetherapy.com](http://www.t-onetherapy.com)

# T-ONE MEDI PRO

T-One Medi Pro is the professional electrostimulator with 4 channels for training and rehabilitation

 **Device designed and certified in Italy**

## SPECIFICATIONS

- ✓ 4 independent channels (up to 8 pluggable electrode)
- ✓ Constant current generator
- ✓ Compensated biphasic square wave
- ✓ Monophasic wave (only for iontophoresis programs)
- ✓ Two-user function (2+2)
- ✓ Backlit color display
- ✓ Operation with 4.8 Volt rechargeable battery pack
- ✓ Adjustable intensity from 1 to 120 mA
- ✓ Wave frequency 0.2 to 1200 Hz
- ✓ Pulse width from 20  $\mu$ s to 300 ms
- ✓ Time from 1 to 90 minutes
- ✓ **Includes treatments for denervation and iontophoresis**



**INCLUDED WITH THE EQUIPMENT**



## PROGRAMS

- |                      |                                 |
|----------------------|---------------------------------|
| ✓ 16 REHA programs   | ✓ 9 URO programs                |
| ✓ 16 TENS programs   | ✓ 22 double user programs (2+2) |
| ✓ 20 NEMS programs   | ✓ 12 TENS/NEMS free memories    |
| ✓ 18 BEAUTY programs |                                 |

For more information about the device and its programs: [www.t-onetherapy.com](http://www.t-onetherapy.com)

## T-ONE REHAB

T-One Rehab is the 4-channel electrostimulator designed for rehabilitation specialists and clinics.

 **Device designed and certified in Italy**

### SPECIFICATIONS

- ✓ 4 independent channels (up to 8 pluggable electrode)
- ✓ Constant current generator
- ✓ Compensated biphasic square wave
- ✓ Monophasic wave (only for iontophoresis programs)
- ✓ Two-user function (2+2)
- ✓ Backlit color display
- ✓ Operation with 4.8 Volt rechargeable battery pack
- ✓ Adjustable intensity from 1 to 120 mA
- ✓ Wave frequency 0.2 to 1200 Hz
- ✓ Pulse width from 20  $\mu$ s to 300 ms
- ✓ Time from 1 to 90 minutes
- ✓ **Includes treatments for denervation and iontophoresis**



**INCLUDED WITH THE EQUIPMENT**



### PROGRAMS

- |                              |                              |
|------------------------------|------------------------------|
| ✓ 9 URO Programs             | ✓ 20 NEMS Programs           |
| ✓ 12 TENS/NEMS free memories | ✓ 22 Two-user programs (2+2) |
| ✓ 16 TENS Programs           | ✓ 29 REHA Programs           |
| ✓ 18 BEAUTY Programs         |                              |

For more information about the device and its programs: [www.t-onetherapy.com](http://www.t-onetherapy.com)

## ELJ41X414/ELJ48X484

4 adhesive electrodes  
pre-gelled with  
2mm jack.  
Sizes: 41x41mm  
48x48mm



**ELJ41X414**  
(41x41mm)

**ELJ48X484**  
(48x48mm)



## ELJ40X804/ELJ50X904

4 adhesive electrodes  
pre-gelled with  
2mm jack.  
Sizes: 40x80mm  
50x90mm



**ELJ40X804**  
(40x80mm)

**ELJ50X904**  
(50x90mm)



## OVAL03

2 adhesive electrodes  
pre-gelled 130x76 mm with  
2mm jack connection.



## ELJ50X1004

2 adhesive electrodes  
pre-gelled 50x100 mm with  
2 jack connections 2mm.



## ELJ30R8

8 adhesive face electrodes  
pre-gelled 30mm with 2mm  
jack connection.



## ELB46X474

4 adhesive electrodes  
pre-gelled 46x47mm with  
button connection.




## ELB53X964

4 adhesive electrodes  
pre-gelled 53x96mm with  
button connection.



## ELB32R8

8 adhesive face electrodes  
pre-gelled 30mm with  
button connection.



## ELADCL8F

4 cable adapters with  
2mm male plug and  
electrode/band with button  
connection



## FP235-245/FC250-265

4 conductive elastic bands for thighs 2x65cm+2x50cm and for  
calves 2x35cm+2x45cm adjustable with button connection.



**FP235-245**  
(2x35cm+2x45cm)

**FC250-265**  
(2x65cm+2x50cm)



**FI80X60**

Blue elastic band  
8x60cm for iontophoresis.



8 019781 471260

**ELB60X801/ELB80X1201**

Deerskin iontophoresis  
electrode with button  
connection.  
Sizes: 60x80mm  
80x120mm



**ELB60X801**  
**(60x80mm)**



8 019781 471277

**ELB80X1201**  
**(80x120mm)**



8 019781 471284

**ELG50X502/ELG60X852****ELG80X1202**

2 rubber electrodes with  
2mm jack connection.  
Sizes: 50x50mm  
60x85mm/80x120mm



**ELG50X502**  
**(50x50mm)**



8 019781 472380

**ELG60X852**  
**(60x85mm)**



8 019781 472397

**ELG80X1202**  
**(80x120mm)**



8 019781 472403

**SPX50X502/SPX60X852****SPX80X1202**

2 x rubber electrode sponges  
with 2mm jack connection.  
Sizes: 50x50mm  
60x85mm/80x120mm



**SPX50X502**  
**(50x50mm)**



8 019781 472410

**SPX60X852**  
**(60x85mm)**



8 019781 472427

**SPX80X1202**  
**(80x120mm)**



8 019781 472434

**KIT-IONOS**

2 electrodes with 2mm jack  
connection, 2 sponges and  
1 elastic band.



8 019781 471604

**IPP-PAD**

Iontophoresis electrode  
specific for IPP (Induratio Penis  
Plastica) treatment.



8 019781 471628

**SA2S/BEAC**

Beac anal probe  
with 2mm jack connection.



8 019781 471369

**SV2S/BEAC**

Beac vaginal probe  
with 2mm jack connection.



8 019781 471376

**SA2S/WHITE**

White anal probe  
with 2mm jack connection.



8 019781 472717

**SV2S/WHITE**

White vaginal probe  
with 2mm jack connection.



8 019781 472724

**GE000250**

Gel 260ml for ultrasound therapy and electrotherapy.



8 019781 470454

**HANDY TRODE**

Handle/Bipolar electrode for face/décolté/breast applications with 2mm jack connection.



8 019781 002020

**PEN-EL**

Point search pen with 2mm jack connection.



8 019781 002013

**CV025318**

Blue and yellow cable for Mio-Care series.



8 019781 471130

**TWINJ2X2**

Splitter cables 2 yellow and 2 blue.



8 019781 471147

**MIO-BATTERY NEW**

Mio-Care and Mag1000 series rechargeable battery (2012)



8 019781 472731

**UNICHAR**

Mio-Care Battery Tens/Beauty/Fitness and Mag1000 (2012)



8 019781 471482

**UNICHAR NEW**

Battery (6.8Vdc) series T-One, Mag1000 (2018), Mio-Vein, Mio-Care Pro, Mio-Ionotens and I-Tech Physio.



8 019781 472885

**T-ONE BATTERY**

Battery series T-One 2016, 4.8V Ni-Mh 2300mA.



8 019781 002037

**CJ025317 / CJ025314 / CJ025312 / CJ025315**

Blue, yellow, red, white cable T-One series with 2mm jack connection. Sold individually



**CJ025317**  
(blue cable)



8 019781 472298



**CJ025314**  
(yellow cable)



8 019781 472274



**CJ025312**  
(red cable)



8 019781 472267



**CJ025315**  
(white cable)



8 019781 472281

# ULTRASOUND THERAPY

Ultrasound is sound waves that travel through the tissues of our body transferring energy. Their passage causes an increase in temperature that promotes microcirculation, reducing pain. If subcutaneous fluids are present, ultrasound promotes its correct elimination.

In addition, the passage of the sound wave allows better penetration into the body of active ingredients useful for therapy.

## THERAPEUTIC BENEFITS

- ✓ Reduction of pain in a few sessions
- ✓ Relaxation of contracted muscles
- ✓ Biostimulation and trophic action

### IDEAL FOR


- All diseases of the musculoskeletal system due to the analgesic effect • Sciatica • Neuritis
- Periarticular calcifications • Ulcers • Duplay's disease • Dupuytren's disease • Tendinitis
- Organized haematoma • Scar tissues • Muscle contractures

Ultrasound therapy is also used in the **aesthetic field** for the treatment of imperfections, the reactivation of local circulation and the absorption of active ingredients.



# POWERSONIC

Powersonic is the new ultrasound therapy device that works both with battery and mains power.

 **Device designed and certified in Italy**

## SPECIFICATIONS

- ✓ 13 treatments
- ✓ Maximum power: 6.4W (with connected power supply) and 3.2W (with battery)
- ✓ 1 MHz ultrasound
- ✓ Maximum surface intensity: 1.6W/cm<sup>2</sup> (with connected power supply) and 0.8W/cm<sup>2</sup> (with battery)
- ✓ 3 selectable intensity levels
- ✓ Collimated ultrasonic beam
- ✓ Operation with medical and battery power supply
- ✓ Suitable for outdoor use (no immersion)

## TREATMENTS

Myofascial pain	Pressure ulcers	Epicondylitis (Golfer's elbow)	Sciatica
Muscle pain	Lumbar pain	Tendinitis	
Carpal tunnel syndrome	Arthrosis	Cubital tunnel syndrome	
Venous ulcers	Epicondylitis	Lumbar stenosis	




**RECOMMENDED  
BY DOCTORS**

**WORKS WIRED AND WIRELESS**

**THE MOST POWERFUL IN ITS CLASS**

## MIO-SONIC

Mio-Sonic is the compact and powerful ultrasonic device, ideal for treating pain, muscle contractures and joint stiffness.

 **Device designed and certified in Italy**

### PROGRAMS

- ✓ 9 settable programs

### SPECIFICATIONS

- ✓ Maximum power: 6.4 W (1.6 W/cm<sup>2</sup>)
- ✓ Ultrasound 1 MHz
- ✓ Collimated ultrasonic beam
- ✓ 3 selectable intensity levels
- ✓ Suitable for outdoor use (no immersion)



# DIAGNOSTICS

The pulse oximeter is an essential tool for measuring the **level of oxygen saturation in the blood** and the heartbeat. Its use is simple and intuitive. It can be used anywhere, even at home

## FOX-220

### SPECIFICATIONS


- ✓ LED Displays
- ✓ Oxygen saturation indication (SpO<sub>2</sub>), heartbeat, pulse bar, perfusion index
- ✓ Auto shutdown when finger is not inserted
- ✓ Low battery indication
- ✓ Fits the finger size of adults and children



## FOX-250

### SPECIFICATIONS


- ✓ Detection and indication of the degree of saturation of blood oxygen and heartbeat
- ✓ Pulse bar indication
- ✓ LED Displays
- ✓ Auto shutdown when finger is not inserted
- ✓ Low battery indication
- ✓ Fits the finger size of adults and children



## FOX-350

### SPECIFICATIONS


- ✓ Detection and indication of the degree of blood oxygen saturation, heartbeat, perfusion index and plethismogram
- ✓ Pulse bar indication
- ✓ OLED display with brightness adjustment
- ✓ Auto shutdown when finger is not inserted
- ✓ Low battery indication
- ✓ Fits the finger size of adults and children



## Digital Infrared Thermometer FI06

### SPECIFICATIONS

- ✓ Fever Alarm Indicator
- ✓ 3-color backlit LCD display
- ✓ Instant Read
- ✓ Accurate measurement
- ✓ Front measurement in degrees C°/F°
- ✓ Available memory



## Surgical masks

### SPECIFICATIONS

- ✓ Facial mask with 3-ply protection
- ✓ Excellent wearability




## FFP2 masks

### SPECIFICATIONS

- ✓ Adaptable nasal band
- ✓ Anti-smog and dustproof
- ✓ Individually wrapped masks

### REATHERM

Resistive and capacitive diathermy.



### LA500

Low power laser therapy up to 500 mW.



### LA8000/LA10000

High power laser therapy up to 8 W/10 W.



### I-TECH UT2

Professional ultrasound therapy 1/3 MHz, double handpiece (5 cm<sup>2</sup> - 1 cm<sup>2</sup>).



### I-TECH UE

Professional combined therapy (ultrasound + electrotherapy).



### I-TECH PHYSIO 4

4-channel electrotherapy device.



### I-TECH PHYSIO EMG

Electromyography 2 channels + Electrotherapy 4 channels.



### POWER Q-6000 PLUS

Professional pressotherapy with 6 compression chambers.



### TROLL-M/L

3 tier trolley with 4 casters and 1 power outlet.



RETAILER



[www.itechmedicaldivision.com](http://www.itechmedicaldivision.com)




**IACER Srl** - Via S. Pertini, 24/A | 30030 Martellago (VE) - Italy  
Ph.: (+39) 041/5401356 | Tax Code and VAT No.: IT00185480274 | [iacer@iacer.it](mailto:iacer@iacer.it)



Official Supplier



Get the most out of your therapy!  
**Download our new app and start your therapeutic journey.**



Last update: December 2021 | Information content reserved for medical personnel.

BR08-02