

T-ONE MEDI PRO

T-One Medi Pro is the professional 4-channel electrotherapy device for training and rehabilitation.

KIT

- ✓ T-One Medi Pro
- ✓ 4 connection cables (8 connectable electrodes)
- ✓ Self-adhesive square electrodes
- ✓ Self-adhesive rectangular electrodes
- ✓ Self-adhesive round electrodes
- ✓ Iontophoresis kit
- ✓ Ni-MH rechargeable battery
- ✓ Battery charger
- ✓ User manual
- ✓ Electrodes position manual
- ✓ Carrying bag



MEDICAL DEVICE CLASS IIa

T-ONE MEDI PRO is a medical device CE0068

TECHNICAL FEATURES

- ✓ Four independent channels unit (up to 8 connectable electrodes)
- ✓ Constant current generator (can guarantee a constant current even when skin condition changes, the device automatically stops in the event of accidental electrodes disconnection)
- ✓ Compensated biphasic square wave: current quantity from positive pole to negative one is constant to avoid polarization effect
- ✓ Monophasic wave generator (only for iontophoresis programs);
- ✓ 2+2 users programs
- ✓ Backlight graphic colour display
- ✓ User feedback: the device recognizes patient connection
- ✓ Soft-touch easy to use keyboard
- ✓ On - line guide on display for fast instructions;
- ✓ 4,8 V rechargeable battery pack;
- ✓ Personalized programs for man and woman;
- ✓ Adjustable intensity from 1 to 120 mA
- ✓ Frequency 0,2 ÷ 1200 Hz
- ✓ Width impulse 20 μ s ÷ 300 ms
- ✓ Time 1÷90 minutes

DIMENSIONS

150 x 90 x 35 mm

T-ONE MEDI PRO

PROGRAMS

16 REHA PROGRAMS	Plantar stimulation	Firming up	2+2 USERS
Iontophoresis L	Epicondylitis	Toning up	Tens rapid
Iontophoresis M	Epitroclea	Definition	Tens endorphinic
Iontophoresis H	Periarthritis	Modelling	Firming up
MENS Micro-current	Tens Burst	Microlifting	Toning up
Haematomas	Micro-current	Lipolysis abdomen	Definition
Oedemas	20 NEMS PROGRAMS	Lipolysis thighs	Modelling
TENS sequential	Warming up	Lipolysis glutes e hips	Microlifting
Tens Burst	Resistance	Lipolysis arms	Muscle relaxant
Atrophy prevention	Resistant strength	Tissue elasticity	Massage
Atrophy	Basic strength	Capillarisation	Deep capillarisation
Denervated 50ms	Fast strength	Heaviness in legs	Resistant strength
Denervated 100ms	Explosive strength	Microlifting face with handpiece	Basic strength
Denervated 200ms	Deep capillarisation	Microlifting seno with handpiece	Fast strength
Denervato 200ms	Muscle recovery	Microlifting décolleté with handpiece	Explosive strength
Denervated 250ms	Agonist / Antagonist	Face wrinkles with handpiece	Hypertrophy
Denervated 300ms	Sequential tonic contractions	Forehead wrinkles with handpiece	Face microlifting with handpiece
16 TENS PROGRAMS	Sequential phasic contractions	Collagen	Breast microlifting with handpiece
TENS rapid	Muscle relaxant	9 URO PROGRAMS	Cleavage microlifting with handpiece
TENS Endorphinic	Deep massage	Stress incontinence	Facial wrinkles with handpiece
TENS maximum values	EMS	Stress incontinence	Forehead wrinkles with handpiece
Anti-inflammatory	30' Resistance	Stress incontinence	FREE MEMORIES
Neck pain	30' Resistant strength	Urge incontinence	5 Tens mem
Backache/Sciatica	20' Basic strength	Urge incontinence	5 NEMS mem
Sprains / Bruises	20' Fast strength	Urge incontinence	NEMS Alternated
Vascularisation	20' Explosive strength	Mixed incontinence	NEMS Alternated
Muscle relaxant	20' Hypertrophy	Mixed incontinence	Battery test
Hand/wrist pain	18 BEAUTY PROGRAMS	Mixed incontinence	