

F2S700

The 2-solenoids therapeutic belt F2S700 is a CE medical device compatible with the high intensity/low frequency pulsed electromagnetic field (PEMF) devices.

Compatible with:

- ✓ MAG700

CLASS I MEDICAL DEVICE, non-sterile and without measuring functions



TECHNICAL FEATURES

- ✓ Applicator size: 38x9 cm + elastic of 60 cm
- ✓ 2 solenoids
- ✓ 1 connection cable

KIT

- ✓ F2S700
- ✓ User Manual