

F3S2000

The 3-solenoids therapeutic belt F3S2000 is a CE medical device compatible with the high intensity/low frequency pulsed electromagnetic field (PEMF) devices.

Compatible with:

- ✓ MAG2000
- ✓ MAG2000 Plus
- ✓ MAG2000 Premium
- ✓ LaMagneto
- ✓ LaMagneto Pro
- ✓ LaMagneto X

CLASS I MEDICAL DEVICE, non-sterile and without measuring functions



TECHNICAL FEATURES

- ✓ Applicator size: 38x9 cm + elastic of 60 cm
- ✓ 3 solenoids
- ✓ 1 connection cable

KIT

- ✓ F3S2000
- ✓ User Manual