

# MIO-CARE PRO

Mio-Care Pro is the ideal electrotherapy device for pain therapy and muscle stimulation.

## TECHNICAL FEATURES

- ✓ Combined electrotherapy unit with 2 independent channels
- ✓ Compensated biphasic square wave: current quantity from positive pole to negative one is constant to avoid hazardous thermal effect of polarization
- ✓ TENS - NEMS - BEAUTY programs
- ✓ 13 free memories TENS - NEMS
- ✓ Digital easy to use keyboard
- ✓ Remaining therapy time and programs instructions shown on graphic display
- ✓ Internal rechargeable battery pack (removable)
- ✓ Backlight display
- ✓ Maximum intensity 200 mApp



### MEDICAL DEVICE CLASS IIa

MIO-CARE PRO is a medical device CE0068

## KIT

- ✓ Mio-Care Pro
- ✓ TENS and NEMS electrotherapy with graphic display
- ✓ 2 connection cables with splitters (8 electrodes connection)
- ✓ Pre-gelled adhesive electrodes
- ✓ Internal Ni-Mh rechargeable battery pack with new LL technology (long lasting)
- ✓ Battery charger
- ✓ Belt clip
- ✓ User manual and electrodes positions manual
- ✓ Carriage bag

## DIMENSIONS

114 X 68 X 26 mm

# MIO-CARE PRO

## PROGRAMS

### 14 TENS PROGRAMS

Conventional tens (fast)

Endorphinic tens (delayed)

Tens at maximum values

Anti-inflammatory

Neck pain/cervicogenic headache

Backache/sciatic pain

Sprains/bruises

Vascularization

Muscle relaxant

Hand and wrist pain

Plantar stimulation

Epicondylitis

Epitroclea

Periarthritis

### 21 NEMS PROGRAMS

Warming up

Resistance - upper limbs and trunk

Resistance - lower limbs

Resistant strength - upper limbs and trunk

Resistant strength - lower limbs

Basic strength - upper limbs and trunk

Basic strength - lower limbs

Fast strength - upper limbs and trunk

Fast strength - lower limbs

Explosive strength - upper limbs and trunk

Explosive strength - lower limbs

Deep capillarization

Muscle recovery

Agonist-antagonist

Sequential tonic contractions  
upper limbs and trunk

Sequential tonic contractions  
lower limbs

Sequential phasic contractions  
upper limbs and trunk

Sequential phasic contractions  
lower limbs

Muscle relaxant

Deep massage

EMS

### 15 BEAUTY PROGRAMS

Warming up

Firming up - upper limbs and trunk

Firming up - lower limbs

Toning up - upper limbs and trunk

Toning up - lower limbs

Definition - upper limbs and trunk

Definition - lower limbs

Modelling

Microlifting

Lipolysis - abdomen

Lipolysis - thighs

Lipolysis - glutei and hips

Lipolysis - arms

Tissue elasticity

Capillarization

Heaviness in legs

### 13 FREE MEMORIES

5 free memories TENS

7 free memories NEMS

1 battery TEST