

OSTEOMAT



OSTEOMAT is the new mattress for low frequency magnetotherapy treatments. This mattress has been realized to meet two important requirements: the long-term therapy and the need of total body treatment.

The configuration of Osteomat stimulates the body zones (shoulders, column, pelvis) that are more affected by bone mass reduction. This is typical of osteoporosis pathology.

Class I medical device non sterile and without measurement functions

MEDICAL DEVICE I non sterile and without measurement functions

KIT

- ✓ Osteomat
- ✓ User manual

TECHNICAL FEATURES

- ✓ 4 active sectors with 8 solenoids in total
- ✓ 2 connection cables
- ✓ Suitable for MAG2000, MAG2000 PLUS and LaMagneto family models

DIMENSIONS

160 x 65 x 2 cm