

WPB-202

Professional device for paraffin baths suitable for hands, feet and the whole body. Paraffin wax is particularly indicated for arthrosis, rheumatoid arthritis and posttraumatic edemas treatments thanks to its effective thermal and mechanical action.

Paraffin is very efficient in preparing the manoeuvres for restoring and treating limbs following immobilization in plaster and posttraumatic stiffness.

The use of paraffin for beauty treatments is more recent. Paraffin has proven to have not only a relief effect on painful joints but also a smoothing effect on the skin. Actually, paraffin moistures skin while opening the pores and facilitating the absorption of active or regenerative ingredients contained by various lotions or masks.



MEDICAL DEVICE CLASS IIa

WPB-202 is a medical device CE0068

BIOLOGICAL AND THERAPEUTIC EFFECTS:

- ✓ **THERMAL EFFECT:** paraffin used at temperatures higher than 50 °C transfers heat to underlying tissues, promoting an exogenous thermal therapeutic action.
- ✓ **MECHANICAL EFFECT:** paraffin applied on skin converts from liquid to solid state. During this transformation, its volume reduces exerting a positive effect by pressing the tissues, particularly an anti-swelling effect.

KIT

- ✓ WPB-202
- ✓ 453 gr x 4 paraffin wax packages
- ✓ 1 hand cover
- ✓ 1 foot cover
- ✓ Single-use Vinyl Pack
- ✓ Safety plate

TECHNICAL FEATURES

Power supply	220-240 VAC, 50-60 Hz
Power consumption	240 W
Temperature	4 levels
Temperature selection	Knob
Indicator	LED

DIMENSIONS

410 x 245 x 250 mmt