

FOX-200

Finger pulse is a very important device used to measure blood oxygen saturation level and cardiac pulsation.

Very easy and intuitive to use, with finger pulse you can monitor your physiological functions everywhere: hospitals, medical centres, first aids, dentists or home.



KIT

- ✓ Device with full colour graphic display
- ✓ 2 AAA size batteries
- ✓ Safe hang lace
- ✓ User and service manual

MEDICAL DEVICE CLASS IIa

Fox-200 is a medical device CE0482

TECHNICAL FEATURES

- ✓ Blood oxygen saturation level measure
- ✓ Cardiac rate measure
- ✓ LED display
- ✓ SpO₂ indication, pulse rate indication and pulsation bar
- ✓ Automatic power off after 16 second without finger insert
- ✓ Low battery indication
- ✓ Paediatric (ages four and up) and adult (weight 15-110 kg) use
- ✓ Weight 50 gr (included batteries)

SPO₂

MEASUREMENT RANGE

35-100%

RESOLUTION

1%

MEASURE ACCURACY

from 70 to 100% → ±2%
under 70% → not defined

CARDIAC PULSATION

MEASUREMENT RANGE

25-250 bpm

RESOLUTION

1 bpm

MEASURE ACCURACY

±2bpm

DIMENSIONS AND WEIGHT:

DEVICE

2x37x32mm

PACKAGING BOX

78x60x36mm

INDIVIDUAL PACKAGING WEIGHT

66 gr