

# FOX-300

Finger pulse is a very important device used to measure blood oxygen saturation level and cardiac pulsation.

Very easy and intuitive to use, with finger pulse you can monitor your physiological functions everywhere: hospitals, medical centres, first aids, dentists or home.



## KIT

- ✓ Device with full colour graphic display
- ✓ 2 AAA size batteries
- ✓ Safe hang lace
- ✓ User and service manual

**MEDICAL DEVICE CLASS IIa**  
Fox-300 is a medical device CE0123

## TECHNICAL FEATURES

- ✓ Blood oxygen saturation level measure
- ✓ Cardiac rate measure
- ✓ Two colours OLED display
- ✓ 10 levels brightness regulation
- ✓ SpO<sub>2</sub> indication, pulse rate indication, pulsation bar and plethysmogram
- ✓ 6 different reading positions
- ✓ Automatic power off
- ✓ Low battery indication
- ✓ Paediatric and adult use
- ✓ Weight 33 gr (excluding batteries)

## SPO<sub>2</sub>

### MEASUREMENT RANGE

70-99%

### RESOLUTION

1%

### MEASURE ACCURACY

from 80 to 99% → ±3%  
from 70 to 80% → ±3%

## CARDIAC PULSATION

### MEASUREMENT RANGE

30-235 bpm

### RESOLUTION

1 bpm

### MEASURE ACCURACY

from 30 to 99bpm → +-2bpm  
from 100 to 235bpm → +-2%

## DIMENSIONS AND WEIGHT

### DEVICE

58x32x37mm

### PACKAGING BOX

83x51x61mm

### INDIVIDUAL PACKAGING WEIGHT

87 gr

# FOX-300

**MANUFACTURER**

Beijing Choice Electronic Technology Co., Ltd.  
Room 4104, No. A12 Yuquan Road, Haidian District, 100143 Beijing, PEOPLE'S REPUBLIC OF CHINA

---

**EC REPRESENTATIVE:**

Shanghai International Holding Corporation GmbH (Europe), Eiffestraße 80, 20537 Hamburg, Germany

---

**DISTRIBUTED BY**

I.A.C.E.R Srl Via S.Pertini 24/a - 30030 Martellago (VE) Italy