

MIO-CARE PRO

Mio-Care Pro is the ideal electrotherapy device for pain therapy and muscle stimulation.

TECHNICAL FEATURES

- ✓ Combined electrotherapy unit with 2 independent channels
- ✓ Compensated biphasic square wave: current quantity from positive pole to negative one is constant to avoid hazardous thermal effect of polarization
- ✓ TENS - NEMS - BEAUTY programs
- ✓ 13 free memories TENS - NEMS
- ✓ Digital easy to use keyboard
- ✓ Remaining therapy time and programs instructions shown on graphic display
- ✓ Internal rechargeable battery pack (removable)
- ✓ Backlight display
- ✓ Maximum intensity 200 mApp



MEDICAL DEVICE CLASS IIa

MIO-CARE PRO is a medical device CE0068

KIT

- ✓ Mio-Care Pro
- ✓ TENS and NEMS electrotherapy with graphic display
- ✓ 2 connection cables with splitters (8 electrodes connection)
- ✓ Pre-gelled adhesive electrodes
- ✓ Internal Ni-Mh rechargeable battery pack with new LL technology (long lasting)
- ✓ Battery charger
- ✓ Belt clip
- ✓ User manual and electrodes positions manual
- ✓ Carriage bag

DIMENSIONS

114 X 68 X 26 mm

MIO-CARE PRO

PROGRAMS

14 TENS PROGRAMS

Conventional tens (fast)
 Endorphinic tens (delayed)
 Tens at maximum values
 Anti-inflammatory
 Neck pain/cervicogenic headache
 Backache/sciatic pain
 Sprains/bruises
 Vascularization
 Muscle relaxant
 Hand and wrist pain
 Plantar stimulation
 Epicondylitis
 Epitroclea
 Periarthritis

13 FREE MEMORIES

5 free memories TENS
 7 free memories NEMS
 1 battery TEST

21 NEMS PROGRAMS

Warming up
 Resistance - upper limbs and trunk
 Resistance - lower limbs
 Resistant strength - upper limbs and trunk
 Resistant strength - lower limbs
 Basic strength - upper limbs and trunk
 Basic strength - lower limbs
 Fast strength - upper limbs and trunk
 Fast strength - lower limbs
 Explosive strength - upper limbs and trunk
 Explosive strength - lower limbs
 Deep capillarization
 Muscle recovery
 Agonist-antagonist
 Sequential tonic contractions
 upper limbs and trunk
 Sequential tonic contractions
 lower limbs
 Sequential phasic contractions
 upper limbs and trunk
 Sequential phasic contractions
 lower limbs
 Muscle relaxant
 Deep massage
 EMS

15 BEAUTY PROGRAMS

Warming up
 Firming up - upper limbs and trunk
 Firming up - lower limbs
 Toning up - upper limbs and trunk
 Toning up - lower limbs
 Definition - upper limbs and trunk
 Definition - lower limbs
 Modelling
 Microlifting
 Lipolysis - abdomen
 Lipolysis - thighs
 Lipolysis - glutei and hips
 Lipolysis - arms
 Tissue elasticity
 Capillarization
 Heaviness in legs