

# POWER Q-1000 PREMIUM

**Power Q-1000 Premium** is the new pressotherapy system with 4 compression chambers. It is indicated both for home use and professional use thanks to its high performance features:


- ✓ Adjustable pressure for each compression chamber
- ✓ Parameters memorization
- ✓ Adjustable air speed level (5,7,9,11,13,15,17,19,21 sec)
- ✓ Adjustable interval time level (5,10,15,20,25,30 sec)
- ✓ 6 setting modes of air chambers (mode A1, mode A2, mode B, mode C, mode D, mode E(A1+C)).

**Power Q-1000 Premium** promotes limb blood flow by the repeated inflation/deflation massage of air cuff. It increases muscle tissues pressure that is indispensable for edema treatment and lower/upper limbs pathologies.



## MEDICAL DEVICE CLASS IIa

Power Q1000 is a medical device CE1023



## TECHNICAL FEATURES

Rated voltage	220-240 V 50/60 Hz
Power consumption	53W
Adjustable time	5÷99 min.
Adjustable pressure	0÷240 mmHg ±20%
Graphic display	LCD
Programs	6
Air speed levels	9 levels
Time interval	6 intervals
Weight	2 Kg

## MAIN APPLICATIONS

- ✓ LYMPHEDEMA / HEMATOMA TREATMENT
- ✓ PASSIVE COMPRESSION FOR VENOUS DISORDERS
- ✓ RHEUMATOID ARTHRITIS
- ✓ ARM DRAINAGE AFTER MASTECTOMY
- ✓ TREATING HEMATOMAS AFTER SURGERIES
- ✓ ANTI-CELLULITE TREATMENT IN THE INITIAL STAGES
- ✓ POST OPERATION PROPHYLAXIS IN VENOUS THROMBOSIS SURGER

# POWER Q-1000 PREMIUM

## KIT

Power Q-1000 Premium can be supplied in 5 different versions:

### POWER Q-1000 Premium TOT

compressive limb therapy system comprising machine body, double legging, 2 derivations connection tube, 2 pressure plantars, armband, abdominal band, 1 derivation connection tube, special connector for simultaneous therapy on both legs and abdominal band.

### POWER Q-1000 Premium LEG2-ABD

compressive limb therapy system comprising machine body, double legging, 2 derivations connection tube, 2 pressure plantars, abdominal band, 1 derivation connection tube, special connector for simultaneous therapy on both legs and abdominal band.

### POWER Q-1000 Premium LEG2

compressive limb therapy system comprising machine body, double legging, 2 derivations connection tube, 2 pressure plantars.

### POWER Q-1000 Premium LEG1

compressive limb therapy system comprising machine body, single legging, 1 derivation connection tube and 1 plantar.

### POWER Q-1000 Premium ARM1

compressive limb therapy system comprising machine body, armband, 1 derivation connection tube.

Extensions for wide legs and special connector for simultaneous therapy on both legs and abdominal band are available as accessories..

## ACCESSORIES SIZE

### Legging

- ✓ Length around 76 cm (from upper side to heel)
- ✓ Thigh circumference around 65 cm
- ✓ Ankle circumference around 38 cm

### Armband

- ✓ Length around 86 cm
- ✓ Biceps circumference around 50 cm
- ✓ Wrist circumference around 36 cm

### Abdominal band

- ✓ Length around 136 cm
- ✓ Height around 46 cm
- ✓ Waistline around 55-110 cm (small) or 80-125 cm (large)



# POWER Q-1000 PREMIUM

**MANUFACTURER**

WONJIN MULSAN Co., Ltd.  
89 Geomdan-ro , Seo-gu , Incheon, Korea

---

**EC REPRESENTATIVE:**

CMC Medical Devices & Drugs S.L. - C/Horacio Lengo N 18, CP 29006, Málaga, Spain

---

**DISTRIBUTED BY**

I.A.C.E.R Srl Via S.Pertini 24/a - 30030 Martellago (VE) Italy