

# POWERSONIC

The POWERSONIC ultrasound-therapy device is ideal for the treatment of musculoskeletal system pathologies, in order to ensure an analgesic and anti-inflammatory effect in problems of various kinds.

The user of the device can be both the patient himself and a professional operator.



## KIT

- ✓ Portable device
- ✓ Medical power supply
- ✓ Ni-Mh AAA850mAh 4.8 V battery
- ✓ Ultrasound gel
- ✓ User manual
- ✓ Carriage and storage bag

**MEDICAL DEVICE CLASS IIa**

Powersonic is a medical device CE0068

## TECHNICAL FEATURES

- ✓ Professional device for ultrasound therapy
- ✓ 3 programs
- ✓ 1MHz ultrasound, maximum output power 6.4W (with power supply) and 3.2W (with battery)
- ✓ Maximum intensity 1.6W/cm<sup>2</sup> (with power supply) and 0.8W/cm<sup>2</sup> (with battery)
- ✓ Collimated ultrasound beam for an optimal penetration beam
- ✓ 3 ultrasound beam intensity levels: L (low), M (medium) e H (high)
- ✓ Medical power supply and battery mode
- ✓ Direct contact device (not suitable for water immersion)

## TREATMENTS

Myofascial pain	Low back pain	Ulnar nerve entrapment
Muscle pain	Arthrosis	Lumbar stenosis
Carpal tunnel	Epicondylitis	Sciatica
Venous ulcers	Golfer's elbow	
Pressure ulcers	Tendinopathy	

## DIMENSIONS

204x63x58mm