



GIMA

BLUETOOTH INFRARED AND EAR THERMOMETER

Code:	25615
Category:	Ear and forehead thermometers
Unit of sale:	1 pc.
Minimum order:	1
Type:	Medical device
Class:	II A
NSIS:	2306429
CND:	V0301010202



Description: **BLUETOOTH INFRARED FOREHEAD AND EAR THERMOMETER**

Multifunctional infrared thermometer used to read body (ear and forehead), room and object temperature.

- clinically accurate, reliable and quick
- professional and home use
- patent technology to separately measure adults and children
- four colour LED backlight
- 10 memory records
- switchable between °C and °F
- audio and visual fever warning
- magnetic probe cover to automatically switch between forehead and ear mode

Bluetooth connection to mobile devices, downloading free Vihealth App.

Supplied with pouch, 2 AAA batteries, multilingual manual: GB, FR, IT, ES, DE.

Technical Specifications:	Measurement type: 3 in 1 (forehead, ear, object, liquid)
	Measurement range/accuracy: forehead & ear 32.0-42.9°C ± 0.2°C / object 0 -100°C / ± 1.0°C
	Auto power off time: 1 min
	Measuring time: 1 second
	Temperature mode: °C/°F switchable
	Memory: 10 groups
	Display: LCD, four-colour LED display (white, green, orange, red)
	Resolution: 0.1 °C/°F
	Battery: 3V DC 2xAAA batteries, can be used more than 3000 times
	Size - weight: 143x35x41 mm - 80g (with batteries)
	Classification - protection: Type BF - IP22