

# SIZE GUIDE

PRESSOTHERAPY

**I-TECH**  
MEDICAL DIVISION

## SIZE GUIDE

### LEGGINS

Leggings with 4 overlapping compression chambers

**THIGH CIRCUMFERENCE**

**ANKLE CIRCUMFERENCE**

**LEG LENGTH**

(measurement taken from inner thigh to heel)

SIZE **M**



Up to 64 cm

Up to 50 cm

Up to 82 cm

SIZE **L**




Up to 70 cm

Up to 50 cm

Up to 98 cm

### LEGGING EXTENDERS



### TNT PROTECTION

Pair of hygienic TNT leg protections



### PRESSURE ORTHOTIC



## SIZE GUIDE

### ARM CUFF

4-compression chamber arm cuff

**ARMPIT CIRCUMFERENCE**

**WRIST CIRCUMFERENCE**

SIZE **M**



Up to 47 cm

Up to 32 cm

SIZE **L**



Up to 56 cm

Up to 40 cm

## NON-MEDICAL KITS

### BEAUTY KIT

- n.1 abdominal band
- n. 1 4-chamber 4-ending tube
- n.1 splitter



### TOTAL KIT

- n.1 abdominal band
- n. 1 arm cuff
- n. 1 4-chamber 4-ending tube
- n.1 splitter



The Total Kit is available with arm cuff in sizes M and L.

## SIZE GUIDE

### ABDOMINAL BAND

4-chamber abdominal band

WAIST CIRCUMFERENCE

SIZE (M)



Up to 116 cm

## TUBES AND CONNECTORS

### 4-CHAMBER 4-ENDING TUBE

Suitable for leg, armband and abdominal band.



### CONNECTOR

4-tube connector for armband, legging or abdominal band.



### SPLITTER

Adapter for simultaneous therapy.  
(2 leggings and abdominal/arm band)



The photos of the accessories are purely illustrative, the product may show variations in some aesthetic details.