



# GIMA

## GIMAFIT BODY FAT SCALE with GimaAPP and Bluetooth - black

**Code:** 27094  
**Category:** Body fat analyzer scales  
**Unit of sale:** 1 pc.  
**Minimum order:** 1



**EAN13:** 8023279270945

**Description:** GIMAFIT BODY SCALE with Gima App and Bluetooth 5.0

Multifunction scale can be connected to GIMA App for detection of weight and body fat values.

Communication protocols and SDK (Software development key) available for Telemedicine projects.

Equipped with 4 precision G sensors and platform made of glass and stainless steel.

Multilingual manual in: GB, FR, IT, ES, PT, DE, PL, GR.

**Technical Specifications:**  
Capacity: 180 kg / 396 lb - minimum 0.15 kg  
LED display: 80x40 mm white light  
Sensitivity: 200 g / 0.4 lb (>50 kg <100 kg)  
Batteries 4xAAA (included)  
Size 300x300x23 mm  
Switchable kg/lb