

# LAMAGNETO X

LaMagneto X is the ideal PEMF device for night-time and/or long-lasting treatments. It is particularly indicated for end users suffering from osteoporosis or pathologies and pains spread throughout the whole body.

Medical device certified in accordance with Regulation (EU) 2017/745 (MDR).

## TECHNICAL FEATURES

- ✓ Two independent channels
- ✓ Colour display
- ✓ Magnetic field intensity up to 200 gauss per channel
- ✓ Frequency from 1 to 120 Hz
- ✓ 63 programs available
- ✓ Double user
- ✓ List of the last 10 programs used
- ✓ "Lock" function
- ✓ Counter of therapies performed and total hours of use



**EMDN CODE:** Z120606  
**MEDICAL DEVICE CLASS** IIa  
**BASIC UDI-DI:** 8019781PEMFLFDEVP2  
**UDI-DI LaMagneto X:** 08019781202031



## PROGRAMS

01 Osteoporosis	22 Cartilage lesion	43 Night 15 Hz
02 Arthrosis	23 Ligament lesion	44 Night 50 Hz
03 Arthritis	24 Bone oedema	45 Night 75 Hz
04 Cervical arthrosis	25 Tendinitis	46 Autoscan
05 Articular pain	26 Chondropathy	47 Treatment 1 Hz
06 Cervicalgias	27 Antinflammatory	48 Treatment 4 Hz
07 Chronic pain	28 Whiplash Syndrome	49 Treatment 5 Hz
08 Fractures	29 Healing	50 Treatment 8 Hz
09 Epicondylitis	30 Cutaneous ulcers	51 Treatment 10 Hz
10 Epitrocleitis	31 Discopathy	52 Treatment 15 Hz
11 Pseudoarthrosis	32 Myalgia	53 Treatment 16 Hz
12 Lumbalgy	33 Neuropathy	54 Treatment 20 Hz
13 Lumbar pain	34 Muscle strain	55 Treatment 30 Hz
14 Shoulder arthrosis	35 Muscle cramp	56 Treatment 40 Hz
15 Knee arthrosis	36 Rizarthrosis	57 Treatment 50 Hz
16 Scapulo Humeral Periarthritis	37 Impingement Syndrome	58 Treatment 60 Hz
17 Coxarthrosis	38 Carpal Tunnel	59 Treatment 70 Hz
18 Muscular atrophy	39 Titanium prosthesis	60 Treatment 75 Hz
19 Muscular treatment	40 Rotator cuff	61 Treatment 80 Hz
20 Osteonecrosis	41 Tarsal Tunnel	62 Treatment 90 Hz
21 Algodystrophy	42 Night 1 Hz	63 Treatment 100 Hz

# LAMAGNETO X

## KIT

- ✓ LaMagneto X
- ✓ A 3 solenoids belt applicator
- ✓ A professional applicator with a pair of solenoids
- ✓ Medical power supply
- ✓ User manual
- ✓ Magnet for verifying therapy operation
- ✓ Carrying bag
- ✓ Three non-woven fabric strips 15x150 cm
- ✓ Available in combination with the Osteomat mattress for night or long-lasting body treatments



**DIMENSIONS** 180 x 110 x 50 mm