

LAMAGNETO

LaMagneto is the latest model of PEMF therapy that allows you to reduce pain and inflammation in the comfort of your own home.

It comes with 46 preset therapeutic programs, updated according to the most recent clinical trials, and 2 independent channels to do different treatments on two parts of the body at the same time, even in pairs.

TECHNICAL FEATURES

- ✓ Two independent channels
- ✓ Colour display
- ✓ Magnetic field intensity up to 150 gauss per channel
- ✓ Frequency from 1 to 100 Hz
- ✓ 46 programs available
- ✓ Double user
- ✓ List of the last 10 programs used



EMDN CODE: Z120606

MEDICAL DEVICE CLASS IIa

LaMagneto is a medical device CE 1936

UDI-DI LaMagneto: 08019781101068

MD

PROGRAMS

01 Osteoporosis	17 Coxarthrosis	33 Neuropathy
02 Arthrosis	18 Muscular atrophy	34 Muscle strain
03 Arthritis	19 Muscular treatment	35 Muscular cramp
04 Cervical arthrosis	20 Osteonecrosis	36 Rizarthrosis
05 Articular pain	21 Algodystrophy	37 Impingement syndrome
06 Cervicalgias	22 Cartilage lesion	38 Carpal Tunnel
07 Chronic pain	23 Ligament lesion	39 Titanium prosthesis
08 Fractures	24 Bone oedema	40 Rotator cuff
09 Epicondylitis	25 Tendinitis	41 Tarsal Tunnel
10 Epitrocleitis	26 Chondropathy	42 Night 1 Hz
11 Pseudoarthrosis	27 Antinflammatory	43 Night 15 Hz
12 Lumbalgia	28 Whiplash syndrome	44 Night 50 Hz
13 Lumbar pain	29 Healing	45 Night 75 Hz
14 Shoulder arthrosis	30 Cutaneous ulcers	46 Autoscan
15 Knee arthrosis	31 Discopathy	
16 Scapulo Humeral Periarthritis	32 Myalgia	

LAMAGNETO

KIT

- ✓ LaMagneto
- ✓ Two 3 solenoids belt applicator
- ✓ Medical power supply
- ✓ User manual
- ✓ Magnet for verifying therapy operation
- ✓ Carrying bag
- ✓ Three non-woven fabric strips 15x150 cm

DIMENSIONS 180 x 110 x 50 mm

