



# GIMA

## RI-CHAMPION SMART PRO+ BLOOD PRESSURE MONITOR WITH BLUETOOTH- 1735-BT - adult cuff

**Code:** 49951  
**Category:** Riester digital sphygmomanometers  
**Unit of sale:** 1 pc.  
**Minimum order:** 1  
**Type:** Medical device  
**Class:** II A  
**NSIS:** 2668630  
**CND:** Z1203020501  
**EAN13:** 4045396218895



**Description:** RI-CHAMPION® SMART PRO+ BLOOD PRESSURE MONITOR WITH BLUETOOTH - 1735 BT - adult cuff

A clinically validated blood pressure monitor, designed for fast and reliable readings. Advanced performance with irregular heartbeat detection and averaging capabilities: analysis of 3 successive blood pressure measurements.

- Bluetooth® connectivity
- fast and quiet operation
- 200 measurements on one set of 4 AA batteries
- large, easy to read backlit screen
- meets EN1060-4, NIBP clinical investigation

**Technical Specifications:** Heart Rate Range: 40-199 beats/minute  
Systolic Measurement Range: 60-255 mmHg  
Diastolic Measurement Range: 30-195 mmHg  
Pulse Rate Measurement Range: 40-199 beats / minute  $\pm 4\%$  of reading  
Maximum Inflation Pressure: 280 mmHg  $\pm 3$  mmHg or  $\pm 2\%$  of reading  
Power Source/Supply: 4x 1.5 V AA alkaline batteries / DC + 6 V / 1 A  
Battery Life: 200 Readings  
Memory: maximum 400 memory records  
Size - weight: 141x121xh 72 mm - 350g without battery