

REACARE

ReaCare is the capacitive-resistive diathermy device with treatments for rehabilitation and sports, which can be customised thanks to free memories.

It is ideal for treating joint pain, knee arthrosis and oedema, but also for helping athletes recover after an injury.



TECHNICAL FEATURES

- ✓ Colour display
- ✓ 24 REHAB programs
- ✓ 10 SPORT programs
- ✓ CUSTOM program
- ✓ 10 Free memories
- ✓ Power 160 W, frequency 400÷600 kHz modulation up to 100 Hz

MEDICAL DEVICE IIb

ReaCare is a medical device CE0068

KIT

- ✓ ReaCare device
- ✓ Resistive probe
- ✓ Capacitive probe
- ✓ Capacitive electrodes 60 mm and 35 mm Ø
- ✓ Resistive electrodes 60 mm and 35 mm Ø
- ✓ Return plate
- ✓ High conductive cream

REACARE

PROGRAMS

REHABILITATION

- 01 Phantom limb
- 02 Knee arthrosis
- 03 Arthrosynovitis
- 04 Bursitis
- 05 Neck pain
- 06 Whiplash syndrome
- 07 Chondropathy
- 08 Muscle contracture
- 09 Contusion
- 10 Muscular sprain
- 11 Sprain
- 12 Muscle distraction
- 13 Joint pain

- 14 Oedema
- 15 Epicondylitis
- 16 Epitrochleitis
- 17 Anterior cruciate ligament (ACL)
- 18 Lymphoedema
- 19 Low back pain
- 20 Meniscopathy (meniscal synovitis)
- 21 Muscle tear
- 22 Acute tendinitis
- 23 Chronic tendinitis
- 24 Rotator cuff tendinopathy

SPORT

- 01 Muscle contracture
- 02 Contusion
- 03 Muscle recovery
- 04 Sprain
- 05 Hematoma
- 06 Tendon injury
- 07 Superficial massage
- 08 Deep massage
- 09 Muscle sprain
- 10 Muscle tear