

REATHERM

Reatherm is our latest capacitive-resistive diathermy model, with neuromodulation effect. It is, indeed, the perfect ally both for therapists and for patients. Thanks to its 400 frequencies (300-700 kHz) and to its modulation system (1-100Hz), healthcare professionals can perfectly adapt any treatment according to the patients' needs, obtaining excellent results both for acute pathologies and chronic ones. On the other hand, patients relieve from pain faster, feeling better starting from the first session.



TECHNICAL FEATURES

- ✓ Colour display
- ✓ 35 REHAB programs
- ✓ 5 BEAUTY programs
- ✓ 10 SPORT programs
- ✓ CUSTOM program
- ✓ 10 Free memories
- ✓ Power 200 W, frequency 300÷700 kHz, modulation up to 100 Hz.

MEDICAL DEVICE IIb

Reatherm is a medical device CE0068

KIT

- ✓ Reatherm device
- ✓ Resistive probe
- ✓ Capacitive probe
- ✓ Capacitive electrodes 60 mm and 35 mm Ø
- ✓ Resistive electrodes 60 mm and 35 mm Ø
- ✓ Bipolar probe with electrode 50 mm Ø
- ✓ Return plate
- ✓ High conductive cream

REATHERM

PROGRAMS

REHAB

- 01 Phantom limb
- 02 Knee arthrosis
- 03 Arthrosynovitis
- 04 Bursitis
- 05 Neck pain
- 06 Whiplash syndrome
- 07 Chondropathy
- 08 Muscle contracture
- 09 Contusion
- 10 Muscular sprain
- 11 TMJ dysfunctions (Jaw articulation dysfunction)
- 12 Sprain
- 13 Muscle distraction
- 14 Joint pain
- 15 Oedema
- 16 Enthesitis
- 17 Epicondylitis
- 18 Epitrochleitis

- 19 Slipped disc
- 20 Cicatricial outcome
- 21 Anterior cruciate ligament (ACL)
- 22 Lymphoedema
- 23 Low back pain
- 24 Sciatica
- 25 Meniscopathy (meniscal synovitis)
- 26 Diabetic foot
- 27 Piriformis syndrome
- 28 Muscle tear
- 29 Acute tendinitis
- 30 Chronic tendinitis
- 31 Rotator cuff tendinopathy
- 32 Ulcers
- 33 Effusion (joint)
- 34 Arthroplasty (recovery)
- 35 Fracture

SPORT

- 01 Muscle contracture
- 02 Contusion
- 03 Muscle recovery
- 04 Sprain
- 05 Hematoma
- 06 Tendon injury
- 07 Superficial massage
- 08 Deep massage
- 09 Muscle sprain
- 10 Muscle tear

BEAUTY

- 01 Cellulitis
- 02 Collagen
- 03 Edema
- 04 Draining massage
- 05 Microlifting