

# MIO-CARE PRO

Mio-Care Pro is the ideal electrotherapy device for pain therapy and muscle stimulation.

## TECHNICAL FEATURES

- ✓ Combined electrotherapy unit with 2 independent channels
- ✓ Compensated biphasic square wave: current quantity from positive pole to negative one is constant to avoid hazardous thermal effect of polarization
- ✓ TENS - NEMS - BEAUTY programs
- ✓ 13 free memories TENS - NEMS
- ✓ Digital easy to use keyboard
- ✓ Remaining therapy time and programs instructions shown on graphic display
- ✓ Internal rechargeable battery pack (removable)
- ✓ Backlight display
- ✓ Maximum intensity 200 mApp



**MEDICAL DEVICE CLASS IIa**

MIO-CARE PRO is a medical device CE0068

## KIT

- ✓ Mio-Care Pro
- ✓ TENS and NEMS electrotherapy with graphic display
- ✓ 2 connection cables with splitters (8 electrodes connection)
- ✓ Pre-gelled adhesive electrodes
- ✓ Internal Ni-Mh rechargeable battery pack with new LL technology (long lasting)
- ✓ Battery charger
- ✓ Belt clip
- ✓ User manual and electrodes positions manual
- ✓ Carriage bag

## DIMENSIONS

114 X 68 X 26 mm

# MIO-CARE PRO

## PROGRAMS

### 14 TENS PROGRAMS

Conventional tens (fast)  
 Endorphinic tens (delayed)  
 Tens at maximum values  
 Anti-inflammatory  
 Neck pain/cervicogenic headache  
 Backache/sciatic pain  
 Sprains/bruises  
 Vascularization  
 Muscle relaxant  
 Hand and wrist pain  
 Plantar stimulation  
 Epicondylitis  
 Epitroclea  
 Periarthritis

### 21 NEMS PROGRAMS

Warming up  
 Resistance - upper limbs and trunk  
 Resistance - lower limbs  
 Resistant strength - upper limbs and trunk  
 Resistant strength - lower limbs  
 Basic strength - upper limbs and trunk  
 Basic strength - lower limbs  
 Fast strength - upper limbs and trunk  
 Fast strength - lower limbs  
 Explosive strength - upper limbs and trunk  
 Explosive strength - lower limbs  
 Deep capillarization  
 Muscle recovery  
 Agonist-antagonist  
 Sequential tonic contractions  
 upper limbs and trunk  
 Sequential tonic contractions  
 lower limbs  
 Sequential phasic contractions  
 upper limbs and trunk  
 Sequential phasic contractions  
 lower limbs  
 Muscle relaxant  
 Deep massage  
 EMS

### 15 BEAUTY PROGRAMS

Warming up  
 Firming up - upper limbs and trunk  
 Firming up - lower limbs  
 Toning up - upper limbs and trunk  
 Toning up - lower limbs  
 Definition - upper limbs and trunk  
 Definition - lower limbs  
 Modelling  
 Microlifting  
 Lipolysis - abdomen  
 Lipolysis - thighs  
 Lipolysis - glutei and hips  
 Lipolysis - arms  
 Tissue elasticity  
 Capillarization  
 Heaviness in legs

### 13 FREE MEMORIES

5 free memories TENS  
 7 free memories NEMS  
 1 battery TEST