

T-ONE MEDI PRO

T-One Medi Pro is the professional 4-channel electrotherapy device for training and rehabilitation.



UMDNS 13762

MEDICAL DEVICE CLASS IIa

TECHNICAL FEATURES

- Four independent channels unit (up to 8 connectable electrodes)
- Constant current generator (can guarantee a constant current even when skin condition changes, the device automatically stops in the event of accidental electrodes disconnection)
- Compensated biphasic square wave: current quantity from positive pole to negative one is constant to avoid polarization effect
- Monophasic wave generator (only for iontophoresis programs)
- 2+2 user programs
- Backlight graphic colour display
- User feedback: the device recognizes patient connection
- Soft-touch easy to use keyboard
- On-line guide on display for fast instructions
- 4,8V rechargeable battery pack
- Personalized programs for man and woman
- Adjustable intensity from 1 to 120 mA
- Frequency 0,2 ÷ 1200 Hz
- Width impulse 20µs ÷ 300ms
- Time 1 ÷ 90 minutes
- Dimensions: 150x90x35 mm

T-ONE MEDI PRO

KIT

- T-One Medi Pro device
- 4 connection cables (8 connectable electrodes)
- Self-adhesive square electrodes
- Self-adhesive rectangular electrodes
- Iontophoresis kit
- Ni-MH rechargeable battery
- Battery charger
- User manual
- Electrodes position manual
- Carrying bag



PROGRAMS

16 REHA PROGRAMS

- 01 Ionophoresis L (low)
- 02 Ionophoresis M (medium)
- 03 Ionophoresis H (high)
- 04 MENS Micro-current
- 05 Haematomas
- 06 Oedemas
- 07 TENS sequential
- 08 TENS Burst
- 09 Atrophy prevention
- 10 Atrophy
- 11 Denervated 50ms
- 12 Denervated 100ms
- 13 Denervated 150ms
- 14 Denervated 200ms
- 15 Denervated 250ms
- 16 Denervated 300ms

16 TENS PROGRAMS

- 01 TENS rapid
- 02 TENS Endorphinic
- 03 TENS maximum values
- 04 Anti-inflammatory
- 05 Neck pain
- 06 Backache/Sciatica
- 07 Sprains / Bruises
- 08 Vascularisation
- 09 Muscle relaxant
- 10 Hand/wrist pain
- 11 Plantar stimulation
- 12 Epicondylitis
- 13 Epitroclea
- 14 Periarthritis
- 15 Tens Burst
- 16 Micro-current

20 NEMS PROGRAMS

- 01 Warming up
- 02 Resistance
- 03 Resistant strength
- 04 Basic strength
- 05 Fast strength
- 06 Explosive strength
- 07 Deep capillarisation
- 08 Muscle recovery
- 09 Agonist / Antagonist
- 10 Sequential tonic contractions
- 11 Sequential phasic contractions
- 12 Muscle relaxant
- 13 Deep massage
- 14 EMS
- 15 30' Resistance
- 16 30' Resistance strength
- 17 20' Basic strength
- 18 20' Fast strength
- 19 20' Explosive strength
- 20 20' Hypertrophy

T-ONE MEDI PRO

PROGRAMS

18 BEAUTY PROGRAMS

01	Firming up
02	Toning up
03	Definition
04	Modelling
05	Microlifting
06	Lipolysis abdomen
07	Lipolysis Thighs
08	Lipolysis glutes and hips
09	Lipolysis arms
10	Tissue elasticity
11	Capillarisation
12	Heaviness in legs
13	Microlifting face with handpiece
14	Microlifting breast with handpiece
15	Microlifting décolleté with handpiece
16	Face wrinkles with handpiece
17	Forehead wrinkles with handpiece
18	Collagen

9 URO PROGRAMS

01	Stress incontinence
02	Stress incontinence
03	Stress incontinence
04	Urge incontinence
05	Urge incontinence
06	Urge incontinence
07	Mixed incontinence
08	Mixed incontinence
09	Mixed incontinence

FREE MEMORIES

M1-M5	TENS free memories
M6-M10	NEMS free memories
M11-M12	NEMS alternated free memories
M13	Battery test

22 2+2 USER PROGRAMS

01	Tens rapid
02	Tens endorphinic
03	Firming up
04	Toning up
05	Definition
06	Modelling
07	Microlifting
08	Muscle relaxant
09	Massage
10	Deep capillarisation
11	Anti-inflammatory
12	Resistance
13	Resistant strength
14	Basic strength
15	Fast strength
16	Explosive strength
17	Hypertrophy
18	Face microlifting with handpiece
19	Breast microlifting with handpiece
20	Cleavage microlifting with handpiece
21	Facial wrinkles with handpiece
22	Forehead wrinkles with handpiece