

T-ONE REHAB

T-One Rehab is the highly professional 4-channel electrotherapy device intended for rehabilitation specialists and rehabilitation clinics.



UMDNS 13762

MEDICAL DEVICE CLASS IIa

TECHNICAL FEATURES

- Four independent channels unit (up to 8 connectable electrodes)
- Constant current generator (can guarantee a constant current even when skin condition changes, the device automatically stops in the event of accidental electrodes disconnection)
- Compensated biphasic square wave: current quantity from positive pole to negative one is constant to avoid polarization effect
- Monophasic wave generator (only for iontophoresis programs)
- 2+2 user programs
- Backlight graphic colour display
- User feedback: the device recognizes patient connection
- Soft-touch easy to use keyboard
- On-line guide on display for fast instructions
- 4,8V rechargeable battery pack
- Personalized programs for man and woman
- Adjustable intensity from 1 to 120 mA
- Frequency 0,2 ÷ 1200 Hz
- Width impulse 20µs ÷ 300ms
- Time 1 ÷ 90 minutes
- Dimensions: 150x90x35 mm

T-ONE REHAB

KIT

- T-One Rehab device
- 4 connection cables (8 connectable electrodes)
- Self-adhesive square electrodes
- Self-adhesive rectangular electrodes
- Iontophoresis kit
- Ni-MH rechargeable battery
- Battery charger
- User manual
- Electrodes position manual
- Carrying bag



PROGRAMS

29 REHA PROGRAMS

01	Ionophoresis L (low)
02	Ionophoresis M (medium)
03	Ionophoresis H (high)
04	MENS Micro-current
05	Haematomas
06	Oedemas
07	TENS sequential
08	TENS Burst
09	Atrophy prevention
10	Atrophy
11	Denervated 50ms
12	Denervated 100ms
13	Denervated 150ms
14	Denervated 200ms
15	Denervated 250ms
16	Denervated 300ms
17	Triangular 50ms
18	Triangular 100ms
19	Triangular 150ms
20	Triangular 200ms
21	Triangular 250ms
22	Triangular 300ms
23	Interferential

24	TENS amplitude modulation
25	TENS alternated
26	TENS frequency modulation
27	TENS frequency modulation
28	TENS amplitude modulation
29	Kotz Wave

16 TENS PROGRAMS

01	TENS rapid
02	TENS Endorphinic
03	TENS maximum values
04	Anti-inflammatory
05	Neck pain
06	Backache/Sciatica
07	Sprains / Bruises
08	Vascularisation
09	Muscle relaxant
10	Hand/wrist pain
11	Plantar stimulation
12	Epicondylitis
13	Epitroclea
14	Periarthritis
15	Tens Burst
16	Micro-current

20 NEMS PROGRAMS

01	Warming up
02	Resistance
03	Resistant strength
04	Basic strength
05	Fast strength
06	Explosive strength
07	Deep capillarisation
08	Muscle recovery
09	Agonist / Antagonist
10	Sequential tonic contractions
11	Sequential phasic contractions
12	Muscle relaxant
13	Deep massage
14	EMS
15	30' Resistance
16	30' Resistance strength
17	20' Basic strength
18	20' Fast strength
19	20' Explosive strength
20	20' Hypertrophy

T-ONE REHAB

PROGRAMS

18 BEAUTY PROGRAMS

01	Firming up
02	Toning up
03	Definition
04	Modelling
05	Microlifting
06	Lipolysis abdomen
07	Lipolysis Thighs
08	Lipolysis glutes and hips
09	Lipolysis arms
10	Tissue elasticity
11	Capillarisation
12	Heaviness in legs
13	Microlifting face with handpiece
14	Microlifting breast with handpiece
15	Microlifting décolleté with handpiece
16	Face wrinkles with handpiece
17	Forehead wrinkles with handpiece
18	Collagen

9 URO PROGRAMS

01	Stress incontinence
02	Stress incontinence
03	Stress incontinence
04	Urge incontinence
05	Urge incontinence
06	Urge incontinence
07	Mixed incontinence
08	Mixed incontinence
09	Mixed incontinence

FREE MEMORIES

M1-M5	TENS free memories
M6-M10	NEMS free memories
M11-M12	NEMS alternated free memories
M13	Battery test

22 2+2 USER PROGRAMS

01	Tens rapid
02	Tens endorphinic
03	Firming up
04	Toning up
05	Definition
06	Modelling
07	Microlifting
08	Muscle relaxant
09	Massage
10	Deep capillarisation
11	Anti-inflammatory
12	Resistance
13	Resistant strength
14	Basic strength
15	Fast strength
16	Explosive strength
17	Hypertrophy
18	Face microlifting with handpiece
19	Breast microlifting with handpiece
20	Cleavage microlifting with handpiece
21	Facial wrinkles with handpiece
22	Forehead wrinkles with handpiece