

# ORTHOMAG

Orthomag is the first device for low frequency magnetotherapy (CEMP) to be used anywhere:

- ✓ With battery (in movement).
- ✓ By car with car charger.
- ✓ At home or office with the power supply adapter

The Orthomag applicator is one of a kind. Extremely light, comfortable to put on and take off and FLEXIBLE, it adapts perfectly to the physiology of your body. No low-frequency magnetotherapy applicator on the market has such a reduced thickness and weight. It folds to fit the shape of your body.

It definitely solves the problems of weight, size and flexibility. Orthomag has two programs (50-75Hz) for the treatment of pathologies related to bone tissue and problems related to inflammatory conditions.



## MEDICAL DEVICE CLASS IIa

ORTHOMAG is a medical device certified CE0197

## TECHNICAL CHARACTERISTICS

- ✓ Universal flexible applicator
- ✓ Lithium rechargeable battery
- ✓ 2 programmes (50 Hz – 75 Hz)
- ✓ Up to 25 Gauss per program



## KIT

- ✓ ORTHOMAG
- ✓ Battery charger (approx.1.5 m cable)
- ✓ User manual
- ✓ Universal flexible applicator (approx.1.5 m cable)
- ✓ Carriage bag
- ✓ Test emissions magnet
- ✓ Two elastic bands (S and L size)
- ✓ Car charger (optional)

# ORTHOMAG

## MANUFACTURER

Shenzhen Dongdixin Technology Co., Ltd.  
Floor 1-2, No.3 Building, Fanshen Xusheng  
Industrial Estate Xilixiaobaimang 518108  
Nanshan District, Shenzhen P.R. China

---

## EC REPRESENTATIVE

Shanghai International Holding Corp. GmbH (Europe)  
Eiffestraße 80, 20537 Hamburg, Germany

---

## DISTRIBUTED BY

I.A.C.E.R Srl Via Enzo Ferrari 2 – 30037 Scorzè (VE) Italy