

MIO-CARE FITNESS

Mio-Care Fitness is the ideal electrostimulator for pain treatment, able to meet the needs of the sportsman and indicated for aesthetic purposes.



TECHNICAL FEATURES

- ✓ Combined electrotherapy unit with 2 independent channels
- ✓ Compensated biphasic square wave: current quantity from positive pole to negative one is constant to avoid hazardous thermal effect of polarization
- ✓ Asynchronous channels working (lipolysis, tonic and phasic contraction, agonist/anatgonist)
- ✓ Digital easy to use keyboard
- ✓ Remaining therapy time and programs instructions shown on graphic display
- ✓ Internal rechargeable battery pack (removable)
- ✓ Backlight display
- ✓ Maximum intensity 100 mApp

MEDICAL DEVICE CLASS IIa

MIO-CARE FITNESS is a medical device CE0197

KIT

- ✓ Mio-Care Fitness
- ✓ TENS and EMS electrotherapy with graphic display
- ✓ 2 connection cables with splitters (8 electrodes connection)
- ✓ Pre-gelled adhesive electrodes
- ✓ Internal Ni-Mh rechargeable battery pack with new LL technology (long lasting)
- ✓ Battery charger
- ✓ Belt clip
- ✓ User manual and electrodes positions manual
- ✓ Carriage bag

DIMENSIONS

114 X 68 X 26 mm

MIO-CARE FITNESS

PROGRAMS

20 MEDICAL PROGRAMS

Conventional tens (fast)
 Endorphinic tens (delayed)
 Tens at maximum values
 Anti-inflammatory
 Neck pain/cervicogenic headache
 Backache/sciatic pain
 Distorsioni/contusioni
 Vascularization
 Muscle relaxant
 Haematomas
 Atrophy prevention
 Atrophy (trophism rehabilitation)
 Hand and wrist pain
 Plantar stimulation
 Epicondylitis
 Epitroclea
 Periarthritis
 Microcurrent
 Stress incontinence (probe not incl.)
 Urgency incontinence (probe not incl.)

15 BEAUTY PROGRAMS

Firming up – upper limbs and trunk
 Firming up – lower limbs
 Toning up – upper limbs and trunk
 Toning up – lower limbs
 Definition – upper limbs and trunk
 Definition – lower limbs
 Modelling
 Microlifting
 Lipolysis – abdomen
 Lipolysis – thighs
 Lipolysis – glutei and hips
 Lipolysis – arms
 Tissue elasticity
 Capillarization
 Heaviness in legs

20 FITNESS PROGRAMS

Warming up
 Resistance – upper limbs and trunk
 Resistance – lower limbs
 Resistant strength – upper limbs and trunk
 Resistant strength – lower limbs
 Basic strength – upper limbs and trunk
 Basic strength – lower limbs
 Fast strength – upper limbs and trunk
 Fast strength – lower limbs
 Explosive strength – upper limbs and trunk
 Explosive strength – lower limbs
 Deep capillarization
 Muscle recovery
 Agonist-antagonist
 Sequential tonic contractions upper limbs and trunk
 Sequential tonic contractions lower limbs
 Sequential phasic contractions upper limbs and trunk
 Sequential phasic contractions lower limbs
 Muscle relaxant
 Deep massage

ALTERNATE STIMULATION ON TWO CHANNELS

MIO-CARE can generate alternate impulses on the two channels to produce sequential contractions or stimulations. This features avoid the use of MIO-CARE FITNESS both for beauty treatments (for example lipolysis) and for fitness treatments to stimulate agonist muscles or in sequential contractions. Its high power can satisfy even the most exigent sport user.

MIO-CARE FITNESS

MANUFACTURER

Shenzhen Dongdixin Technology Co., Ltd.
Floor 1-2, No.3 Building, Fanshen Xusheng
Industrial Estate Xilixiaobaimang 518108
Nanshan District, Shenzhen P.R. China

EC REPRESENTATIVE

Shanghai International Holding Corp. GmbH (Europe)
Eiffestraße 80, 20537 Hamburg, Germany

DISTRIBUTED BY

I.A.C.E.R Srl Via Enzo Ferrari 2 – 30037 Scorzè (VE) Italy