

MAG1000

Mag1000 is the high frequency PEMF device particularly recommended for the treatment of pathologies affecting soft tissues like muscles and tendons.



TECHNICAL FEATURES

- ✓ 20 programs
- ✓ 2 output channels
- ✓ Adjustable intensity (3 levels)
- ✓ Carrier frequency 27 MHz
- ✓ Frequency from 8 Hz to 640 Hz

MEDICAL DEVICE CLASS IIa

MAG1000 is a medical device CE0197.

KIT

- ✓ MAG1000
- ✓ elastic therapeutic belt (FT11000)
- ✓ 1 battery pack
- ✓ Battery charger
- ✓ User manual
- ✓ Carriage bag

DIMENSIONS 153 x 91 x 38 mm

PROGRAMS

01 Arthritis-Arthrosis	08 Sprains	15 Epicondylitis
02 Rheumatisms	09 Dislocations	16 Epytrocleitis
03 Pain- Articular prosthesis	10 Traumas	17 Lumbago
04 Herniated disc	11 Bruises	18 Muscular contractures
05 Muscular pain	12 Myalgia	19 2h Medium application
06 Osteoporosis	13 Tendonitis	20 4h Long application
07 Fractures	14 Hematoma	

MAG1000

MANUFACTURER

Shenzhen Dongdixin Technology Co., Ltd.
Floor 1-2, No.3 Building, Fanshen Xusheng
Industrial Estate Xilixiaobaimang 518108
Nanshan District, Shenzhen P.R. China

EC REPRESENTATIVE

Shanghai International Holding Corp. GmbH (Europe)
Eiffestraße 80, 20537 Hamburg, Germany

DISTRIBUTED BY

I.A.C.E.R Srl Via Enzo Ferrari 2 – 30037 Scorzè (VE) Italy