

# T-ONE MEDI SPORT

T-One Medi Sport is the 4 channels electro-stimulator dedicated to regular athletes who want to improve their performances and do rehabilitation independently.

## KIT

- ✓ T-One Medi Sport
- ✓ 4 connection cables (8 connectable electrodes)
- ✓ Self-adhesive square electrodes
- ✓ Self-adhesive rectangular electrodes
- ✓ Self-adhesive round electrodes
- ✓ Iontophoresis kit
- ✓ Ni-MH rechargeable battery
- ✓ Battery charger
- ✓ User manual
- ✓ Electrodes position manual
- ✓ Carrying bag



## MEDICAL DEVICE CLASS IIa

T-ONE MEDI SPORT is a medical device CE0068

## TECHNICAL FEATURES

- ✓ Four independent channels unit (up to 8 connectable electrodes)
- ✓ Constant current generator (can guarantee a constant current even when skin condition changes, the device automatically stops in the event of accidental electrodes disconnection)
- ✓ Compensated biphasic square wave: current quantity from positive pole to negative one is constant to avoid polarization effect
- ✓ Monophasic wave generator (only for iontophoresis programs);
- ✓ 2+2 users programs
- ✓ Backlight graphic colour display
- ✓ User feedback: the device recognizes patient connection
- ✓ Soft-touch easy to use keyboard
- ✓ On – line guide on display for fast instructions;
- ✓ 4,8 V rechargeable battery pack;
- ✓ Personalized programs for man and woman;
- ✓ Adjustable intensity from 1 to 120 mA
- ✓ Frequency 1 ÷ 1200 Hz
- ✓ Width impulse 20 ÷ 450 µs
- ✓ Time 1÷90 minutes

## DIMENSIONS

150 x 90 x 35 mm

# T-ONE MEDI SPORT

## PROGRAMS

<b>10 REHA PROGRAMS</b>	<b>20 NEMS PROGRAMS</b>	Lipolysis thighs	Explosive strength
Iontophoresis L	Warming up	Lipolysis glutes e hips	Hypertrophy
Iontophoresis M	Resistance	Lipolysis arms	Face microlifting with handpiece
Iontophoresis H	Resistant strength	Tissue elasticity	Breast microlifting with handpiece
MENS Micro-current	Basic strength	Capillarisation	Cleavage microlifting with handpiece
Haematomas	Fast strength	Heaviness in legs	Facial wrinkles with handpiece
Oedemas	Explosive strength	Microlifting face with handpiece	Forehead wrinkles with handpiece
TENS sequential	Deep capillarisation	Microlifting seno with handpiece	<b>FREE MEMORIES</b>
Tens Burst	Muscle recovery	Microlifting décolleté with handpiece	5 Tens mem
Atrophy prevention	Agonist / Antagonist	Face wrinkles with handpiece	5 NEMS mem
Atrophy	Sequential tonic contractions	Forehead wrinkles with handpiece	NEMS Alternated
<b>16 TENS PROGRAMS</b>	Sequential phasic contractions	Collagen	NEMS Alternated
Tens rapid	Muscle relaxant	<b>2+2 USERS</b>	Battery test
Tens Endorphinic	Deep massage	Tens rapid	
Tens maximum values	EMS	Tens endorphinic	
Anti-inflammatory	30' Resistance	Firming up	
Neck pain	30' Resistant strength	Toning up	
Backache/Sciatica	20' Basic strength	Definition	
Sprains / Bruises	20' Fast strength	Modelling	
Vascularisation	20' Explosive strength	Microlifting	
Muscle relaxant	20' Hypertrophy	Muscle relaxant	
Hand/wrist pain	<b>18 BEAUTY PROGRAMS</b>	Massage	
Plantar stimulation	Firming up	Deep capillarisation	
Epicondylitis	Toning up	Anti-inflammatory	
Epitroclea	Definition	Resistance	
Periarthritis	Modelling	Resistant strength	
Tens Burst	Microlifting	Basic strength	
Micro-current	Lipolysis abdome	Fast strength	