

T-ONE REHAB

T-One Rehab is the highly professional 4-channel electrotherapy device intended for rehabilitation specialists and rehabilitation clinics.

KIT

- ✓ T-One Rehab
- ✓ 4 connection cables (8 connectable electrodes)
- ✓ Self-adhesive square electrodes
- ✓ Self-adhesive rectangular electrodes
- ✓ Self-adhesive round electrodes
- ✓ Iontophoresis kit
- ✓ Ni-MH rechargeable battery
- ✓ Battery charger
- ✓ User manual
- ✓ Electrodes position manual
- ✓ Carrying bag



MEDICAL DEVICE CLASS IIa

T-ONE REHAB is a medical device CE0068

TECHNICAL FEATURES

- ✓ Four independent channels unit (up to 8 connectable electrodes)
- ✓ Constant current generator (can guarantee a constant current even when skin condition changes, the device automatically stops in the event of accidental electrodes disconnection)
- ✓ Compensated biphasic square wave: current quantity from positive pole to negative one is constant to avoid polarization effect
- ✓ Monophasic wave generator (only for iontophoresis programs);
- ✓ 2+2 users programs
- ✓ Backlight graphic colour display
- ✓ User feedback: the device recognizes patient connection
- ✓ Soft-touch easy to use keyboard
- ✓ On – line guide on display for fast instructions;
- ✓ 4,8 V rechargeable battery pack;
- ✓ Personalized programs for man and woman;
- ✓ Adjustable intensity from 1 to 120 mA
- ✓ Frequency 0,2 ÷ 1200 Hz
- ✓ Width impulse 20 μ s ÷ 300 ms
- ✓ Time 1÷90 minutes

DIMENSIONS

150 x 90 x 35 mm

T-ONE REHAB

PROGRAMS

29 REHA PROGRAMS	TENS Endorphinic	20' Basic strength	Mixed incontinence
Iontophoresis L	TENS maximum values	20' Fast strength	22 2+2 USERS PROGRAMS
Iontophoresis M	Anti-inflammatory	20' Explosive strength	Tens rapid
Iontophoresis H	Neck pain	20' Hypertrophy	Tens endorphinic
MENS Micro-current	Backache / Sciatica	18 BEAUTY PROGRAMS	Firming up
Haematomas	Sprains / Bruises	Firming up	Toning up
Oedemas	Vascularisation	Toning up	Definition
TENS sequential	Muscle relaxant	Definition	Modelling
Tens Burst	Hand / wrist pain	Modelling	Microlifting
Atrophy prevention	Plantar stimulation	Microlifting	Muscle relaxant
Atrophy	Epitroclea	Lipolysis abdomen	Massage
Denervated 50ms	Epitroclea	Lipolysis thighs	Deep capillarisation
Denervated 100ms	Periarthritis	Lipolysis glutes e hips	Anti-inflammatory
Denervated 150ms	Tens Burst	Lipolysis arms	Resistance
Denervated 200ms	Micro-current	Tissue elasticity	Resistant strength
Denervated 250ms	20 NEMS PROGRAMS	Capillarisation	Basic strength
Denervated 300ms	Warming up	Heaviness in legs	Fast strength
Triangular 50ms	Resistance	Microlifting face with handpiece	Explosive strength
Triangular 100ms	Resistant strength	Microlifting breast with handpiece	Hypertrophy
Triangular 150ms	Basic strength	Microlifting décolleté with handpiece	Face microlifting with handpiece
Triangular 200ms	Fast strength	Face wrinkles with handpiece	Breast microlifting with handpiece
Triangular 250ms	Explosive strength	Forehead wrinkles with handpiece	Cleavage microlifting with handpiece
Triangular 300ms	Deep capillarisation	Collagen	Facial wrinkles with handpiece
Interferential	Muscle recovery	9 URO PROGRAMS	Forehead wrinkles with handpiece
TENS amplitude modulation	Agonist / Antagonist	Stress incontinence	FREE MEMORIES
TENS alternated	Sequential tonic contractions	Stress incontinence	5 Tens mem
TENS frequency modulation	Sequential phasic contractions	Stress incontinence	5 NEMS mem
TENS frequency modulation	Muscle relaxant	Urge incontinence	NEMS Alternated
TENS amplitude modulation	Deep massage	Urge incontinence	NEMS Alternated
Kotz Wave	EMS	Urge incontinence	Battery test
16 TENS PROGRAMS	30' Resistance	Mixed incontinence	
TENS rapid	30' Resistant strength	Mixed incontinence	