

TAP2000

TAP2000 is the 40x40 cm carpet with 2 cables of 3 solenoids each (total 6 solenoids). It is a CE medical device compatible with the high intensity/low frequency pulsed electromagnetic field (PEMF) devices.

Compatible with:

- ✓ LaMagneto
- ✓ LaMagneto Pro
- ✓ LaMagneto X
- ✓ MAG2000
- ✓ MAG2000 Premium
- ✓ MAG2000 Plus

EMDN CODE: Z12069099

CLASS I MEDICAL DEVICE, non-sterile and without measuring functions



TECHNICAL FEATURES

- ✓ Applicator size: 40x40 cm
- ✓ 6 solenoids
- ✓ 2 connection cables

KIT

- ✓ TAP2000
- ✓ User Manual